

添付資料 3

ブナシメジの品種識別マニュアルの
妥当性検査データ(品種識別検査)

ブナシメジ品種識別検査結果票の照合用マーカーフラグメントサイズ一覧表

1. マーカーテーブル

品種番号	品種名	マーカー							
		msHm-1	msHm-2	msHm-3	msHm-4	msHm-5	msHm-6	msHm-7	msHm-8
		A	B	C	D	E	F	G	H
1	ホクト白1号菌	A	B-12/B	C	D	E	F	G	H
2	MH025014	A	B-12/B+8	C-6/C	D+6	E+21	F+4	G-40	H-23
3	MH025016	A-18/A	B-12	C-18/C-6	D-6	E+12	F	G	H-23
4	MH025019	A	B-12/B	C-18/C	D-6	E+15	F+8	G-40/G	H-23
5	MH025022	A	B-12	C-18/C-6	D+6	E+15	F/F+8	G-40/G	H-23
6	MH025066	A-18/A	B-12/B	C-6	D	E/E+15	F/F+4	G-40	H
7	MH025091	A	B-12/B	C-18	D-6	E+15	F	G-40/G	H-23
8	MH025146	A	B-12	C	D-6	E/E+15	F+4	G	H-36
9	MH025217	A-18/A	B-12/B	C	D	E/E+15	F	G	H
10	MH025374	A-18/A	B-12/B	C	D	E	F	G-60	H
11	MH025038	A-12	B-12	C-26/C-8	D+6	E+9/E+24	F+4	G-45	H-23
12	MH025045	A-18	B-12/B-8	C-12/C	D/D+6	E+24	F	G-35	H-28
13	MH025046	A-18/A-6	B-20/B-8	C-6	D-6	E+12	F+4	G-40	H-23/H
14	MH025054	A-6	B-8	C-18	D	E+9	F-4	G-40	H-23
15	MH025070	A-6/A+12	B-8	C-6	D-12	E+9	F	G-40	H-36
16	MH025092	A-18	B-12/B-4	C	D-6/D+6	E+15	F+4/F+8	G-45	H-23
17	MH025097	A-18	B-12	C-18/C-6	D+6	E-4/E	F+4/F+18	G-55	H-34
18	MH025121	A-18/A-6	B-16	C-18	D+12	E-4/E+15	F	G-60	H-20
19	MH025160	A-18/A-6	B-16/B-4	C-18	D+12	E+12	F-8	G-48	H+48
20	MH025339	A-18/A-12	B-12/B-8	C-12	D+18	E+15	F-4/F+4	G-40	H-23
21	MH025012	A-18/A	B-12/B	C-6	D	E/E+15	F/F+4	G-40	H
22	MH025064	A-18/A-12	B-12	C-18/C-6	D	E-4/E	F+4	G-60	H-36
23	MH025286	A	B-12/B	C-18/C-6	D	E/E+24	F/F+4	G-60	H
24	MH025288	A-12/A	B-12	C-18	D-6	E-4/E+15	F+12	G-40	H-23
25	MH025362	A-18/A-12	B-12/B	C-6	D+6	E	F+4	G-60	H+4
26	MH025369	A-18/A-12	B-12/B-8	C-6	D+6	E/E+12	F+4	G-60	H-26/H+4
27	MH025397	A-18	B-16/B	C	D-6	E/E+15	F	G-60/G	H
28	MH025462	A-12	B-12/B-8	C-18/C-6	D+6	E	F/F+4	G-60	H-26
29	MH025463	A-12	B-8/B	C-6	D	E/E+24	F+4	G-60	H+4
30	MH025528	A	B-12/B	C-6	D	E/E+15	F/F+4	G-40	H

2. 検査結果照合用マーカーテーブル

品種番号	品種名	マーカー															
		msHm-1		msHm-2		msHm-3		msHm-4		msHm-5		msHm-6		msHm-7		msHm-8	
		133		67	79	145		155		269		130		132		243	
1	ホクト白1号菌	133		67	79	145		155		269		130		132		243	
2	MH025014	133		67	87	139	145	161		290		134		92		220	
3	MH025016	115	133	67		127	139	149		281		130		132		220	
4	MH025019	133		67	79	127	145	149		284		138		92	132	220	
5	MH025022	133		67		127	139	161		284		130	138	92	132	220	
6	MH025066	115	133	67	79	139		155		269	284	130	134	92		243	
7	MH025091	133		67	79	127		149		284		130		92	132	220	
8	MH025146	133		67		145		149		269	284	134		132		207	
9	MH025217	115	133	67	79	145		155		269	284	130		132		243	
10	MH025374	115	133	67	79	145		155		269		130		72		243	
11	MH025038	121		67		119	137	161		278	293	134		87		220	
12	MH025045	115		67	71	133	145	155	161	293		130		97		215	
13	MH025046	115	127	59	71	139		149		281		134		92		220	243
14	MH025054	127		71		127		155		278		126		92		220	
15	MH025070	127	145	71		139		143		278		130		92		207	
16	MH025092	115		67	75	145		149	161	284		134	138	87		220	
17	MH025097	115		67		127	139	161		265	269	134	148	77		209	
18	MH025121	115	127	63		127		167		265	284	130		72		223	
19	MH025160	115	127	63	75	127		167		281		122		84		291	
20	MH025339	115	121	67	71	133		173		284		126	134	92		220	
21	MH025012	115	133	67	79	139		155		269	284	130	134	92		243	
22	MH025064	115	121	67		127	139	155		265	269	134		72		207	
23	MH025286	133		67	79	127	139	155		269	293	130	134	72		243	
24	MH025288	121	133	67		127		149		265	284	142		92		220	
25	MH025362	115	121	67	79	139		161		269		134		72		247	
26	MH025369	115	121	67	71	139		161		269	281	134		72		217	247
27	MH025397	115		63	79	145		149		269	284	130		72	132	243	
28	MH025462	121		67	71	127	139	161		269		130	134	72		217	
29	MH025463	121		71	79	139		155		269	293	134		72		247	
30	MH025528	133		67	79	139		155		269	284	130	134	92		243	

バナシメジ品種識別検査結果票(B社・許容範囲±1)

サンプル 番号	検査	項目	msHm-1		msHm-2		msHm-3		msHm-4		msHm-5		msHm-6		msHm-7		msHm-8		判定品種	正解品種名	品種名一致	一致 フラグメント数	最大 フラグメント数	一致率
(1)-1	1回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1			
		正解値 フラグ1	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				13	16	81.3%
	2回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1			
		正解値 フラグ1	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				10	16	62.5%
(1)-2	1回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1			
		正解値 フラグ1	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0				9	16	56.3%
	2回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1			
		正解値 フラグ1	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0					16	56.3%
(1)-3	1回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1			
		正解値 フラグ1	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0				11	16	68.8%
	2回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1			
		正解値 フラグ1	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0				11	16	68.8%
(1)-4	1回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1			
		正解値 フラグ1	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				10	16	62.5%
	2回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1			
		正解値 フラグ1	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				10	16	62.5%
(1)-5	1回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1			
		正解値 フラグ1	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0				10	16	62.5%
	2回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1			
		正解値 フラグ1	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0				10	16	62.5%
(1)-6	1回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1			
		正解値 フラグ1	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0				12	16	75.0%
	2回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1			
		正解値 フラグ1	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0				12	16	75.0%
(1)-7	1回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1			
		正解値 フラグ1	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0				11	16	68.8%
	2回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1			
		正解値 フラグ1	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0				11	16	68.8%
(1)-8	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0			
		正解値 フラグ1	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				11	16	68.8%
	2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0			
		正解値 フラグ1	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				11	16	68.8%
(1)-9	1回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1			
		正解値 フラグ1	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0				12	16	75.0%
	2回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1			
		正解値 フラグ1	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0				12	16	75.0%
(1)-10	1回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1			
		正解値 フラグ1	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0				9	16	56.3%
	2回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1			
		正解値 フラグ1	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0				9	16	56.3%
(1)-11	1回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1			
		正解値 フラグ1	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0				10	16	62.5%
	2回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1			
		正解値 フラグ1	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0				10	16	62.5%

(1)－12	1回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1					
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0								
		フラグ1		1	1	1	1	1	1	1	1	1	1	1	1	0	0	1	1					12		16
2回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1						
	正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0									
	フラグ1		1	1	1	1	1	1	1	1	1	1	1	1	0	0	1	1					12		16	75.0%
(1)－13	1回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1					
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0								
		フラグ1		0	0	1	1	1	1	1	1	1	1	1	0	1	1	1					10		16	62.5%
2回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1						
	正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0									
	フラグ1		0	0	1	1	1	1	1	1	1	1	1	0	1	1	1					10		16	62.5%	
(1)－14	1回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1					
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0								
		フラグ1		0	1	1	1	1	1	1	1	1	1	1	0	1	1	1					11		16	68.8%
2回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1						
	正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0									
	フラグ1		0	1	1	1	1	1	1	1	1	1	1	0	1	1	1									
(1)－15	1回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1					
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0								
		フラグ1		1	1	1	1	1	1	1	1	1	1	1	1	0	0	1	1					12		16
2回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1						
	正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0									
	フラグ1		1	1	1	1	1	1	1	1	1	1	1	1	0	0	1	1					12		16	75.0%
(1)－16	1回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1					
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0								
		フラグ1		1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1					12		16
2回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1						
	正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0									
	フラグ1		1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1					12		16	75.0%
(1)－17	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0					
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0								
		フラグ1		0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1					11		16
2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0						
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0									
	フラグ1		0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1					11		16	68.8%
(1)－18	1回目	測定結果	123		66	71	127	140	150	162	270		135		70											

(1)ー24	1回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	0	1				11	16	68.8%
(1)ー25	2回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	0	1				11	16	68.8%
(1)ー26	1回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	0	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1				10	16	62.5%
(1)ー27	2回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	0	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1				10	16	62.5%
(1)ー28	1回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1				12	16	75.0%
(1)ー29	2回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1				12	16	75.0%
(1)ー30	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1				11	16	68.8%
(1)ー31	2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1				11	16	68.8%
(1)ー32	1回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	0	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1				10	16	62.5%
(1)ー33	2回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	0	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1				10	16	62.5%
(1)ー34	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1				11	16	68.8%
(1)ー35	2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1				11	16	68.8%
(1)ー36	1回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	0	1	1				11	16	68.8%
(1)ー37	2回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	0	1	1				11	16	68.8%
(1)ー38	1回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	0	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1				10	16	62.5%
(1)ー39	2回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	0	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1				10	16	62.5%
(1)ー40	1回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー41	2回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%

(1)－36	1回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1				
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1					12	16	75.0%
(1)－37	2回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1				
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1					12	16	75.0%
(1)－38	1回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1				
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0							
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	0	1	1					11	16	68.8%
(1)－39	2回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1				
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0							
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	0	1	1					11	16	68.8%
(1)－40	1回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1				
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	0	1					11	16	68.8%
(1)－41	2回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1				
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	0	1					11	16	68.8%
(1)－42	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0				
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1					11	16	68.8%
(1)－43	2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0				
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1					11	16	68.8%
(1)－44	1回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1				
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1					12	16	75.0%
(1)－45	2回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1				
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1					12	16	75.0%
(1)－46	1回目	測定結果	123		66	71	140		162		270	281	135		70		218	247	MH025286	MH025286	1				
		正解値	121	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247							
		フラグ1	0	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1					10	16	62.5%
(1)－47	2回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1				
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247							
		フラグ1	0	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1					10	16	62.5%
(1)－48	1回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1				
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0							
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1					11	16	68.8%
(1)－49	2回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1				
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0							
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1					11	16	68.8%
(1)－50	1回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1				
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0							
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1					12	16	75.0%
(1)－51	2回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1				
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0							
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1					12	16	75.0%
(1)－52	1回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1				
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0							
		フラグ1	0	1	1	1	1	1	1	0	1	1	0	1	0	1	1	1					9	16	56.3%
(1)－53	2回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1				
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0							
		フラグ1	0	1	1	1	1	1	1	0	1	1	0	1	0	1	1	1					9	16	56.3%
(1)－54	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0				
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1					11	16	68.8%
(1)－55	2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0				
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1					11	16	68.8%

(1)－48	1回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1					
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						11	16	68.8%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1								
(1)－49	1回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1					
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						11	16	68.8%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1								
(1)－49	2回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1					
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						10	16	62.5%
		フラグ1	0	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1								
(1)－50	1回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1					
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						10	16	62.5%
		フラグ1	0	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1								
(1)－50	2回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1					
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1								
(1)－51	1回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1					
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1								
(1)－51	2回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1					
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						10	16	62.5%
		フラグ1	0	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1								
(1)－52	1回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1					
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1								
(1)－52	2回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1					
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1								
(1)－53	1回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1					
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						11	16	68.8%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1								
(1)－53	2回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1					
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						11	16	68.8%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1								
(1)－54	1回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1					
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						11	16	68.8%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	0	1								
(1)－54	2回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1					
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						11	16	68.8%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	0	1								
(1)－55	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0					
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						11	16	68.8%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1								
(1)－55	2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0					
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						11	16	68.8%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1								
(1)－56	1回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1					
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						10	16	62.5%
		フラグ1	0	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1								
(1)－56	2回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1					
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						10	16	62.5%
		フラグ1	0	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1								
(1)－57	1回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1					
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						6	16	37.5%
		フラグ1	1	1	1	0	0	1	1	1	0	0	0	0	1	1	1	1								
(1)－57	2回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1					
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						6	16	37.5%
		フラグ1	1	1	1	0	0	1	1	1	0	0	0	0	1	1	1	1								
(1)－58	1回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1					
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						10	16	62.5%
		フラグ1	0	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1								
(1)－58	2回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1					
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						10	16	62.5%
		フラグ1	0	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1								
(1)－59	1回目	測定結果	123		66	71	127	140	150	162	270		13													

(1)ー60	1回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	0	1	1				11	16	68.8%
	2回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	0	1	1				11	16	68.8%
		totalフラグ計	14	44	60	60	60	60	54	60	60	54	60	6	50	56	60		正解 品種数	100			平均 マーカ ー一致率	67.2%
			58		120		120		114		120		114		56		116							
		一致率	24%		50%		50%		48%		50%		48%		23%		48%							

ブナシメジ品種識別検査結果票(F社・許容範囲±1)

サンプル 番号	検査	項目	msHm-1		msHm-2		msHm-3		msHm-4		msHm-5		msHm-6		msHm-7		msHm-8		判定品種	正解品種名	品種名一致	一致 フラグメント数	最大 フラグメント数	一致率
(1)ー1	1回目	測定結果	117	122	66	71	139	0	161	0	270	281	131	135	70	0	217	247	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				13	16	81.3%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	0	0	1	1	1						
(1)ー2	2回目	測定結果	117	122	66	71	139	0	161	0	270	281	135	0	70	0	217	246	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				14	16	87.5%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1						
(1)ー3	1回目	測定結果	122	146	66	71	127	139	149	162	270	0	131	135	70	0	217	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0				13	16	81.3%
		フラグ1	1	0	1	1	1	1	1	0	1	1	1	1	0	1	1	1						
(1)ー4	2回目	測定結果	122	146	66	71	127	139	149	161	270	0	131	135	70	0	217	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0				13	16	81.3%
		フラグ1	1	0	1	1	1	1	1	0	1	1	1	1	0	1	1	1						
(1)ー5	1回目	測定結果	122	134	66	0	127	145	149	155	266	284	131	0	70	90	217	247	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0				9	16	56.3%
		フラグ1	1	1	1	1	1	0	1	0	1	1	0	1	0	0	0	0						
(1)ー6	2回目	測定結果	122	134	66	0	127	0	149	0	266	284	131	143	70	90	×	0	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0				12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	0	0	0	1	1						
(1)ー7	1回目	測定結果	117	122	66	70	139	0	162	0	270	281	131	135	70	0	217	247	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				13	16	81.3%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	0	0	1	1	1						
(1)ー8	2回目	測定結果	117	122	66	71	139	0	161	0	270	281	135	0	70	0	217	247	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				14	16	87.5%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1						
(1)ー9	1回目	測定結果	117	122	66	79	127	139	162	0	270	0	135	0	70	0	247	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0				13	16	81.3%
		フラグ1	0	1	1	1	1	0	1	1	1	1	1	1	1	1	1	1						
(1)ー10	2回目	測定結果	117	122	66	79	139	0	161	0	270	0	135	0	70	0	246	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0				14	16	87.5%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1	1						
(1)ー11	1回目	測定結果	134	0	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1						
(1)ー12	2回目	測定結果	134	0	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1						
(1)ー13	1回目	測定結果	117	0	×	0	146	0	149	0	270	284	131	0	×	0	243	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0				11	16	68.8%
		フラグ1	0	1	0	0	0	1	1	1	1	1	1	1	0	0	1	1						
(1)ー14	2回目	測定結果	117	0	62	79	146	0	149	0	270	284	131	0	70	0	243	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0				13	16	81.3%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	0	1	1						
(1)ー15	1回目	測定結果	117	134	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025066・MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				14	16	87.5%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1						
(1)ー16	2回目	測定結果	117	134	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025066・MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				14	16	87.5%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1						
(1)ー17	1回目	測定結果	134	146	66	79	127	139	156	0	270	292	131	0	70	0	243	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0				13	16	81.3%
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	0	0	1	1	1						
(1)ー18	2回目	測定結果	134	146	67	79	127	139	155	0	270	292	131	0	70	0	242	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0				13	16	81.3%
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	0	0	1	1	1						
(1)ー19	1回目	測定結果	122	146	66	71	127	139	149	161	270	0	131	135	70	0	217	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0				13	16	81.3%
		フラグ1	1	0	1	1	1	1	1	0	1	1	1	1	0	1	1	1						
(1)ー20	2回目	測定結果	117	122	66	71	127	139	149	161	270	0	131	135	70	0	217	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0				13	16	81.3%
		フラグ1	1	0	1	1	1	1	1	0	1	1	1	1	0	1	1	1						
(1)ー21	1回目	測定結果	117	122	66	0	127	139	156	0	266	270	131	135	70	0	206	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0				13	16	81.3%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	0	0	1	1	1						
(1)ー22	2回目	測定結果	117	122	66	0	127	139	156	0	266	270	135	143	70	0	×	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0				12	16	75.0%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	0	0	1	0	1						

(1)ー12	1回目	測定結果	134	0	66	79	145	0	156	0	270	0	131	0	70	132	243	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1						
2回目	測定結果	134	0	66	79	145	0	155	0	270	0	131	0	×	0	243	0	ホクト白1号菌	ホクト白1号菌	1				
	正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0				15	16	93.8%	
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1							
(1)ー13	1回目	測定結果	116	122	66	79	139	0	161	0	270	0	135	0	70	0	246	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1						
2回目	測定結果	117	122	66	79	139	0	161	0	270	0	135	0	70	0	246	0	MH025362	MH025362	1				
	正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0				14	16	87.5%	
	フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1							
(1)ー14	1回目	測定結果	122	146	71	79	139	0	156	0	270	292	131	135	70	0	247	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0				13	16	81.3%
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	0	1	1	1	1						
2回目	測定結果	122	146	71	79	139	0	155	0	×	0	135	0	70	0	247	0	MH025463	MH025463	1				
	正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0				12	16	75.0%	
	フラグ1	1	0	1	1	1	1	1	1	0	0	1	1	0	1	1	1							
(1)ー15	1回目	測定結果	134	0	66	79	145	0	155	0	270	0	131	0	70	132	243	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1						
2回目	測定結果	134	0	66	79	145	0	155	0	270	0	131	0	70	132	242	0	ホクト白1号菌	ホクト白1号菌	1				
	正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0				15	16	93.8%	
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1							
(1)ー16	1回目	測定結果	134	0	66	79	127	139	155	0	270	292	131	0	70	0	243	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0				14	16	87.5%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	0	0	1	1	1						
2回目	測定結果	134	146	67	79	127	139	156	0	270	292	131	135	70	0	243	0	MH025286	MH025286	1				
	正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0				14	16	87.5%	
	フラグ1	1	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1							
(1)ー17	1回目	測定結果	×	0	66	79	139	0	155	0	270	284	131	135	90	0	243	0	MH025066・MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				13	16	81.3%
		フラグ1	0	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1						
2回目	測定結果	117	134	66	79	139	0	155	0	270	284	131	135	90	0	243	0	MH025066・MH025012	MH025066	0				
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				14	16	87.5%	
	フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1							
(1)ー18	1回目	測定結果	122	146	66	71	127	139	149	161	270	0	131	135	70	0	217	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0				13	16	81.3%
		フラグ1	1	0	1	1	1	1	1	0	1	1	1	1	0	1	1	1						
2回目	測定結果	122	146	66	71	127	139	149	161	270	0	131	135	70	0	217	0	MH025462	MH025462	1				
	正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0				13	16	81.3%	
	フラグ1	1	0	1	1	1	1	1	0	1	1	1	1	0	1	1	1							
(1)ー19	1回目	測定結果	117	134	66	79	139	0	155	0	270	284	131	135	90	0	242	0	MH025066・MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				14	16	87.5%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1						
2回目	測定結果	117	134	66	79	139	0	156	0	270	284	131	135	90	0	242	0	MH025066・MH025012	MH025012	0				
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				14	16	87.5%	
	フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1							
(1)ー20	1回目	測定結果	117	134	66	79	139	0	155	0	270	284	131	135	70	90	246	0	MH025066・MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				10	16	62.5%
		フラグ1	1	0	1	1	1	1	1	0	1	1	0	1	0	1	0	1						
2回目	測定結果	117	134	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025066・MH025012	MH025012	0				
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				10	16	62.5%	
	フラグ1	1	0	1	1	1	1	1	0	1	1	0	1	0	1	0	1							
(1)ー21	1回目	測定結果	117	122	66	71	139	145	156	161	270	281	131	0	70	0	217	246	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				6	16	37.5%
		フラグ1	0	0	1	0	1	0	0	1	1	0	1	0	0	1	0	0						
2回目	測定結果	117	122	66	71	139	0	161	0	270	281	135	0	70	132	217	247	MH025369	MH025369	1				
	正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				7	16	43.8%	
	フラグ1	0	0	1	0	1	1	0	1	1	0	1	0	0	1	0	0							
(1)ー22	1回目	測定結果	134	0	66	79	145	0	156	0	270	0	131	0	70	132	243	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0				12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1						
2回目	測定結果	134	0	66	79	145	0	155	0	270	0	131	0	70	132	243	0	ホクト白1号菌	ホクト白1号菌	1				
	正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0				12	16	75.7	

(1)－36	1回目	測定結果	134	0	66	79	145	0	156	0	270	284	131	135	90	0	242	0	MH025528	MH025528	1				
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0					14	16	87.5%
		フラグ1	1	1	1	1	1	0	1	1	1	1	1	1	0	1	1	1							
2回目	測定結果	134	0	66	79	139	0	156	0	270	284	131	135	90	0	242	0	MH025528	MH025528	1					
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1					15	16	93.8%
(1)－37	1回目	測定結果	117	0	62	79	146	0	149	0	270	284	131	0	70	0	242	0	MH025397	MH025397	1				
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0							
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	0	0	1	1					13	16
2回目	測定結果	116	0	62	79	146	0	149	0	270	284	131	0	70	0	243	0	MH025397	MH025397	1					
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0							
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	0	1	1	1					14	16
(1)－38	1回目	測定結果	117	122	66	0	127	0	155	0	266	284	131	0	90	0	×	0	MH025288	MH025288	1				
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
		フラグ1	1	0	1	1	1	1	1	0	1	1	1	0	1	0	1	0	1					11	16
2回目	測定結果	122	134	66	0	127	0	149	0	266	284	131	143	70	90	243	0	MH025288	MH025288	1					
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	0	0	0	0	1	1					12	16
(1)－39	1回目	測定結果	111	117	66	79	139	145	156	0	270	284	131	135	90	0	242	0	MH025066・MH025012	MH025012	0				
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	0	0	1	1	1	0	1	1	1	1	1	1	0	1	1	1	1					12	16
2回目	測定結果	117	134	66	79	139	0	155	0	270	284	131	135	90	0	243	0	MH025066・MH025012	MH025012	0					
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1	1					14	16
(1)－40	1回目	測定結果	134	0	66	79	139	0	156	0	270	284	135	0	70	90	242	0	MH025528	MH025528	1				
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	1	1	1	1	1	1	1	1	1	1	0	1	0	0	1	1	1					13	16
2回目	測定結果	134	0	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025528	MH025528	1					
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	0	1	1	1					15	16
(1)－41	1回目	測定結果	122	134	66	0	128	0	149	155	266	284	131	0	70	90	206	221	MH025288	MH025288	1				
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
		フラグ1	1	1	1	1	1	1	1	0	1	1	1	0	1	0	0	1	0					11	16
2回目	測定結果	122	134	66	0	127	0	149	155	266	284	143	0	70	90	×	0	MH025288	MH025288	1					
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92										

(1)ー48	1回目	測定結果	122	146	71	79	139	0	156	0	270	292	135	0	70	0	246	0	MH025463	MH025463	1				
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0					14	16	87.5%
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1							
	2回目	測定結果	122	146	71	79	139	0	155	0	270	292	135	0	70	0	246	0	MH025463	MH025463	1				
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0					14	16	87.5%
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1							
(1)ー49	1回目	測定結果	117	122	66	71	139	0	161	0	270	281	135	0	70	0	217	246	MH025369	MH025369	1				
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247					14	16	87.5%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1							
	2回目	測定結果	117	122	66	71	139	0	161	0	270	281	135	0	70	0	217	246	MH025369	MH025369	1				
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247					14	16	87.5%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1							
(1)ー50	1回目	測定結果	134	0	66	79	127	139	155	0	270	292	131	135	70	0	242	0	MH025286	MH025286	1				
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0					15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1							
	2回目	測定結果	134	146	67	79	127	139	156	0	270	292	131	0	70	0	217	242	MH025286	MH025286	1				
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0					12	16	75.0%
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	0	0	1	1	0							
(1)ー51	1回目	測定結果	117	122	66	79	139	0	161	0	270	0	135	0	70	0	246	0	MH025362	MH025362	1				
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0					14	16	87.5%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1							
	2回目	測定結果	116	122	66	79	139	0	161	0	270	0	135	0	70	0	246	0	MH025362	MH025362	1				
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0					15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1							
(1)ー52	1回目	測定結果	122	134	66	79	145	0	156	0	270	0	131	135	70	132	242	0	ホクト白1号菌	ホクト白1号菌	1				
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0					13	16	81.3%
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	0	1	0	1	1							
	2回目	測定結果	122	134	66	79	145	0	156	0	270	0	131	135	70	132	242	246	ホクト白1号菌	ホクト白1号菌	1				
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0					12	16	75.0%
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	0	1	0	1	0							
(1)ー53	1回目	測定結果	122	146	71	79	139	0	156	0	270	292	135	0	70	0	246	0	MH025463	MH025463	1				
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0					14	16	87.5%
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	1	0	1	1	1							
	2回目	測定結果	122	146	71	79	127	139	155	0	270	292	135	0	70	132	246	0	MH025463	MH025463	1				
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0							
		フラグ1	1	0	1	1	1	1	0	1	1	1	1	1	0	0	1	1					12	16	75.0%
(1)ー54	1回目	測定結果	122	134	66	0	127	0	149	168	266	284	131	143	70	90	217	246	MH025288	MH025288	1				
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0					10	16	62.5%
		フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	0	0	0	0							
	2回目	測定結果	122	134	66	0	127	145	149	155	266	284	131	143	70	90	217	0	MH025288	MH025288	1				
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0					10	16	62.5%
		フラグ1	1	1	1	1	1	0	1	0	1	1	1	1	0	0	0	1							
(1)ー55	1回目	測定結果	117	134	66	79	139	0	156	0	270	284	131	135	90	0	217	242	MH025066・MH025012	MH025066	0				
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0					13	16	81.3%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	0							
	2回目	測定結果	117	134	66	79	139	0	155	0	270	284	131	135	90	0	242	0	MH025066・MH025012	MH025066	0				
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1					14	16	87.5%
(1)ー56	1回目	測定結果	117	122	66	79	139	0	161	0	270	0	131	135	70	0	246	0	MH025362	MH025362	1				
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0					13	16	81.3%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1							
	2回目	測定結果	117	122	66	79	139	0	161	0	270	0	135	0	70	0	246	0	MH025362	MH025362	1				
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0					14	16	87.5%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1							
(1)ー57	1回目	測定結果	134	0	66	79	139	0	155	0	270	284	131	135	90	0	217	242	MH025528	MH025528	1				
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0					14	16	87.5%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	0							
	2回目	測定結果	134	0	66	79	139	0	155	0	270	284	131	135	90	0	217	242	MH025528	MH025528	1				
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0					14	16	87.5%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	0							
(1)ー58	1回目	測定結果	117	122	66	0	127	139	156	0	266	270	135	0	70	0	217	0	MH025064	MH025064	1				
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0					13	16	81.3%
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	0	1	0	1							
	2回目	測定結果	117	122	67	0	127	139	156	0															

(1)ー60	1回目	測定結果	117	0	62	79	127	146	149	0	270	284	131	0	70	0	242	0	MH025397	MH025397	1					
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0								
		フラグ1	0	1	1	1	1	0	1	1	1	1	1	1	1	0	0	1	1				12		16	75.0%
	2回目	測定結果	117	0	62	79	146	0	149	0	270	284	131	0	70	0	242	0	MH025397	MH025397	1					
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0								
		フラグ1	0	1	1	1	1	1	1	1	1	1	1	1	1	0	0	1	1				13		16	81.3%
		totalフラグ計	28	35	59	51	56	52	53	48	56	47	54	41	5	44	45	55			正解 品種数	100			平均 マーカ ー致率	79.3%
			63		110		108		101		103		95		49		100									
		一致率	26%		46%		45%		42%		43%		40%		20%		42%									

ブナシメジ品種識別検査結果票(L社・許容範囲±1)

サンプル 番号	検査	項目	msHm-1		msHm-2		msHm-3		msHm-4		msHm-5		msHm-6		msHm-7		msHm-8		判定品種	正解品種名	品種名一致	一致 フラグメント数	最大 フラグメント数	一致率
(1)－1	1回目	測定結果	117	123	69	73	141	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	0	0	0	0	0	1	0	1	0	1	0	1	1	1	0	1				7	16	43.8%
	2回目	測定結果	118	124	69	73	142	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	0	0	0	0	0	1	0	1	0	1	0	1	1	1	0	1				7	16	43.8%
(1)－2	1回目	測定結果	123	0	69	73	128	142	152	165	270	0	136	0	72	0	219	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	0	1	0	0	1	0	0	0	1	1	0	0	1	1	0	1				7	16	43.8%
	2回目	測定結果	124	0	69	73	128	141	152	165	271	0	136	0	72	0	219	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	0	1	0	0	1	0	0	0	1	0	0	0	1	1	0	1				6	16	37.5%
(1)－3	1回目	測定結果	123	135	69	0	128	0	152	0	267	285	145	0	92	0	0	0	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	0	0	0	1	1	1	0	1	0	1	0	1	1	1	0	1				9	16	56.3%
	2回目	測定結果	123	135	69	0	128	0	152	0	267	285	146	0	92	0	0	0	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	0	0	0	1	1	1	0	1	0	1	0	1	1	1	0	1						
(1)－4	1回目	測定結果	118	124	69	73	141	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	0	0	0	0	0	1	0	1	0	1	0	1	1	1	0	1				7	16	43.8%
	2回目	測定結果	118	124	69	73	141	0	165	0	271	282	136	0	71	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	0	0	0	0	0	1	0	1	0	1	0	1	1	1	0	1				7	16	43.8%
(1)－5	1回目	測定結果	117	124	69	81	141	0	165	0	271	0	137	0	72	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	0	0	0	0	0	1	0	1	0	1	0	1	1	1	1	1				8	16	50.0%
	2回目	測定結果	118	124	69	81	141	0	165	0	271	0	136	0	71	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	0	0	0	0	0	1	0	1	0	1	0	1	1	1	1	1				8	16	50.0%
(1)－6	1回目	測定結果	135	0	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	1	0	0	0	1	0	1	0	1	0	0	1	1	1	1				8	16	50.0%
	2回目	測定結果	135	0	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	1	0	0	0	1	0	1	0	1	0	0	1	1	1	1						
(1)－7	1回目	測定結果	118	0	65	81	149	0	152	0	271	285	132	0	72	0	245	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	0	1	0	0	0	1	0	1	0	1	0	1	1	0	0	1				7	16	43.8%
	2回目	測定結果	118	0	65	81	149	0	153	0	271	285	132	0	71	0	244	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	0	1	0	0	0	1	0	1	0	1	0	1	1	0	1	1				8	16	50.0%
(1)－8	1回目	測定結果	118	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	0	0	0	0	1	0	1	0	1	0	0	1	1	1	1				7	16	43.8%
	2回目	測定結果	118	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	0	0	0	0	1	0	1	0	1	0	0	1	1	1	1				7	16	43.8%
(1)－9	1回目	測定結果	135	0	69	81	128	141	159	0	271	293	132	0	72	0	244	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	0	1	0	0	1	0	0	1	0	1	0	0	1	1	1	1				8	16	50.0%
	2回目	測定結果	135	0	69	81	129	141	159	0	271	293	132	0	72	0	244	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	0	1	0	0	0	0	0	1	0	1	0	0	1	1	1	1				7	16	43.8%
(1)－10	1回目	測定結果	123	0	69	73	128	141	152	165	271	0	136	0	72	0	219	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	0	1	0	0	1	0	0	0	1	0	0	0	1	1	0	1				6	16	37.5%
	2回目	測定結果	123	0	69	73	128	142	152	165	270	0	136	0	71	0	219	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	0	1	0	0	1	0	0	0	1	1	0	0	1	1	0	1				7	16	43.8%
(1)－11	1回目	測定結果	117	123	69	0	128	141	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	0	0	0	1	1	0	0	1	0	0	1	1	1	1	0	1				7	16	43.8%
	2回目	測定結果	118	123	69	0	128	142	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	0	0	0	1	1	0	0	1	0	0	1	1	1	1	0	1				7	16	43.8%

(1)ー12	1回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	72	133	244	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	0	1	0	0	0	1	0	1	0	1	0	1	1	0	1	1				8	16	50.0%
(1)ー13	2回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	72	133	244	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	0	1	0	0	0	1	0	1	0	1	0	1	1	0	1	1				8	16	50.0%
(1)ー13	1回目	測定結果	117	124	69	81	141	0	165	0	271	0	136	0	72	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	0	0	0	0	0	1	0	1	0	1	0	1	1	1	1	1				8	16	50.0%
(1)ー13	2回目	測定結果	118	123	69	81	141	0	165	0	271	0	136	0	71	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	0	0	0	0	0	1	0	1	0	1	0	1	1	1	1	1				8	16	50.0%
(1)ー14	1回目	測定結果	123	0	73	81	141	0	159	0	271	293	136	0	72	0	248	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	0	1	0	0	0	1	0	1	0	1	0	1	1	1	1	1				9	16	56.3%
(1)ー14	2回目	測定結果	124	0	73	81	141	0	159	0	271	293	136	0	71	0	248	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	0	1	0	0	0	1	0	1	0	1	0	1	1	1	1	1				9	16	56.3%
(1)ー15	1回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	72	133	244	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	0	1	0	0	0	1	0	1	0	1	0	1	1	0	1	1				8	16	50.0%
(1)ー15	2回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	72	133	244	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	0	1	0	0	0	1	0	1	0	1	0	1	1	0	1	1				8	16	50.0%
(1)ー16	1回目	測定結果	135	0	69	81	128	141	159	0	271	293	132	0	72	0	244	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	0	1	0	0	1	0	0	1	0	0	0	1	1	1	1	1				8	16	50.0%
(1)ー16	2回目	測定結果	135	0	69	81	129	141	159	0	271	293	132	0	72	0	244	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	0	1	0	0	0	0	0	1	0	1	0	0	1	1	1	1				7	16	43.8%
(1)ー17	1回目	測定結果	117	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	0	0	0	0	0	1	0	0	1	0	0	1	1	1	1				7	16	43.8%
(1)ー17	2回目	測定結果	118	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	0	0	0	0	1	0	0	1	0	0	1	1	1	1	1				7	16	43.8%
(1)ー18	1回目	測定結果	123	0	69	72	128	141	152	165	271	0	136	0	72	0	219	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	0	1	0	1	1	0	0	0	0	1	0	0	1	1	0	1				7	16	43.8%
(1)ー18	2回目	測定結果	123	0	69	72	128	141	152	165	271	0	136	0	71	0	219	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	0	1	0	1	1	0	0	0	0	1	0	0	1	1	0	1				7	16	43.8%
(1)ー19	1回目	測定結果	117	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	0	0	0	0	1	0	1	0	1	0	0	1	1	1	1				7	16	43.8%
(1)ー19	2回目	測定結果	118	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	0	0	0	0	1	0	1	0	1	0	0	1	1	1	1				7	16	43.8%
(1)ー20	1回目	測定結果	117	135	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	0	0	0	0	1	0	1	0	1	0	0	1	1	1	1				7	16	43.8%
(1)ー20	2回目	測定結果	118	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	0	0	0	0	1	0	1	0	1	0	0	1	1	1	1				7	16	43.8%
(1)ー21	1回目	測定結果	118	123	69	73	142	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	0	0	0	0	0	1	0	1	0	1	0	1	1	1	0	1				7	16	43.8%
(1)ー21	2回目	測定結果	118	123	69	73	141	0	165	0	271	282	136	0	71	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	0	0	0	0	0	1	0	1	0	1	0	1	1	0	1	1				7	16	43.8%
(1)ー22	1回目	測定結果	135	0	68	81	148	0	159	0	271	0	132	0	72	133	244	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	0	1	1	1	0	0	1	0	1	0	1	1	1	0	1	1				9	16	56.3%
(1)ー22	2回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	71	134	244	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	0	1	0	0	0	1	0	1	0	1	0	1	0	0	1	1				7	16	43.8%
(1)ー23	1回目	測定結果	117	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	0	0	0	0	1	0	1	0	1	0	0	1	1	1	1				7	16	

(1)ー24	1回目	測定結果	123	135	69	0	128	0	152	0	267	285	145	0	92	0	244	0	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	0	0	0	1	1	1	0	1	0	1	0	1	1	1	0	1				9	16	56.3%
(1)ー25	2回目	測定結果	124	135	68	0	128	0	152	0	267	285	145	0	92	0	244	0	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	0	0	0	1	1	1	0	1	0	1	0	1	1	1	0	1				10	16	62.5%
(1)ー26	1回目	測定結果	117	123	69	0	128	141	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	0	0	0	1	1	0	0	1	0	0	0	1	1	1	0	1				7	16	43.8%
(1)ー27	2回目	測定結果	117	123	68	0	128	141	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	0	0	0	1	1	0	0	1	0	0	0	1	1	1	0	1				8	16	50.0%
(1)ー28	1回目	測定結果	135	0	68	81	141	0	159	0	271	285	132	136	92	0	245	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	0	1	1	0	0	1	0	0	1	0	0	1	1	0	1				8	16	50.0%
(1)ー29	2回目	測定結果	135	0	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	0	1	0	0	0	1	0	1	0	0	0	1	1	1	1				8	16	50.0%
(1)ー30	1回目	測定結果	117	135	68	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025066,MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	0	0	1	0	0	1	0	1	0	0	0	1	1	1	1				8	16	50.0%
(1)ー31	2回目	測定結果	117	135	69	81	142	0	159	0	270	285	132	136	92	0	244	0	MH025066,MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	0	0	0	0	1	0	0	1	1	0	0	1	1	1	1				8	16	50.0%
(1)ー32	1回目	測定結果	117	123	68	81	142	0	165	0	271	0	136	0	72	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	0	0	1	0	0	1	0	1	0	1	0	1	1	1	1	1				9	16	56.3%
(1)ー33	2回目	測定結果	118	123	69	81	141	0	165	0	271	0	136	0	71	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	0	0	0	0	0	1	0	1	0	1	0	1	1	1	1	1				8	16	50.0%
(1)ー34	1回目	測定結果	117	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066,MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	0	0	0	0	1	0	0	1	0	0	0	1	1	1	1				7	16	43.8%
(1)ー35	2回目	測定結果	118	135	69	81	142	0	159	0	271	285	132	137	92	0	244	0	MH025066,MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	0	0	0	0	1	0	0	1	0	0	0	1	1	1	1				7	16	43.8%
(1)ー36	1回目	測定結果	118	0	64	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	0	0	1	1	0	0	1	0	1	0	0	1	1	0	1	1				9	16	56.3%
(1)ー37	2回目	測定結果	118	0	65	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	0	0	1	0	0	0	1	0	1	0	0	1	1	0	1	1				8	16	50.0%
(1)ー38	1回目	測定結果	117	123	69	0	128	141	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	0	0	0	0	1	1	0	0	0	0	0	1	1	1	0	1				7	16	43.8%
(1)ー39	2回目	測定結果	117	123	69	0	128	141	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	0	0	0	0	1	1	0	0	0	0	0	1	1	1	0	1				7	16	43.8%
(1)ー40	1回目	測定結果	123	0	72	81	141	0	159	0	271	293	136	0	72	0	248	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	0	0	1	1	0	0	1	0	1	0	0	1	1	1	1	1				10	16	62.5%
(1)ー41	2回目	測定結果	123	0	73	81	141	0	159	0	271	293	136	0	71	0	248	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	0	0	1	0	0	0	1	0	0	1	0	0	1	1	1	1				9	16	56.3%
(1)ー42	1回目	測定結果	117	0	65	81	149	0	152	0	271	285	132	0	72	0	245	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	0	0	1	0	0	0	1	0	1	0	0	1	1	0	0	1				7	16	43.8%
(1)ー43	2回目	測定結果	117	0	65	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	0	0	1	0	0	0	1	0	1	0	0	1	1	0	1	1				8	16	50.0%
(1)ー44	1回目	測定結果	117	123	69	0	128	141	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	0	0	0	0	1	1	0	0	0	0	0	1	1	1	0	1				7	16	43.8%
(1)ー45	2回目	測定結果	117	123	69	0	128	141	159	0	267	271	136	0	71	0	209	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	0	0	0	0	1	1	0	0	0	0	0	1	1	1	0	1				7	16	43.8%
(1)ー46	1回目	測定結果	135	0	68	81	148	0	159	0	271	0	132	0	72	134	244	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	0	0	1	1	0	0	1	0	1	0	0	1	0	0	1	1				8	16	50.0%

(1)ー36	1回目	測定結果	135	0	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1				
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	0	1	0	0	0	1	0	1	0	1	0	0	1	1	1	1				8		16	50.0%
2回目	測定結果	135	0	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1					
	正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0								
	フラグ1	0	1	0	0	0	1	0	1	0	1	0	0	1	1	1	1				8		16	50.0%	
(1)ー37	1回目	測定結果	117	0	64	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1				
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0							
		フラグ1	0	1	1	0	0	1	0	1	0	1	0	1	1	0	1	1				9		16	56.3%
2回目	測定結果	117	0	65	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1					
	正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0								
	フラグ1	0	1	0	0	0	1	0	1	0	1	0	1	1	0	1	1				8		16	50.0%	
(1)ー38	1回目	測定結果	123	135	69	0	128	0	152	0	267	285	145	0	92	0	245	0	MH025288	MH025288	1				
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
		フラグ1	0	0	0	0	1	1	1	0	1	0	1	0	1	1	0	1				9		16	56.3%
2回目	測定結果	123	135	69	0	128	0	152	0	267	285	145	0	92	0	244	0	MH025288	MH025288	1					
	正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0								
	フラグ1	0	0	0	0	1	1	1	0	1	0	1	0	1	1	0	1				9		16	56.3%	
(1)ー39	1回目	測定結果	117	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025012	0				
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	0	0	0	0	0	1	0	1	0	1	0	0	1	1	1	1				7		16	43.8%
2回目	測定結果	117	135	69	81	141	0	159	0	271	285	132	136	92	0	245	0	MH025066.MH025012	MH025012	0					
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0								
	フラグ1	0	0	0	0	0	1	0	1	0	1	0	0	1	1	0	1				6		16	37.5%	
(1)ー40	1回目	測定結果	135	0	69	81	141	0	159	0	271	285	132	136	92	0	245	0	MH025528	MH025528	1				
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	0	1	0	0	0	1	0	1	0	1	0	0	1	1	0	1				7		16	43.8%
2回目	測定結果	135	0	69	81	141	0	159	0	271	285	132	136	92	0	245	0	MH025528	MH025528	1					
	正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0								
	フラグ1	0	1	0	0	0	1	0	1	0	1	0	0	1	1	0	1				7		16	43.8%	
(1)ー41	1回目	測定結果	123	135	69	0	129	0	152	0	267	285	145	0	92	0	244	0	MH025288	MH025288	1				
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
		フラグ1	0	0	0	0	1	0	1	0	1	0	1	0	1	1	0	1				8		16	50.0%
2回目	測定結果	123	135	69	0	129	0	152	0	267	285	145	0	92	0	244	0	MH025288	MH025288	1					
	正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0								
	フラグ1	0	0	0	0	1	0	1	0																

(1)ー48	1回目	測定結果	123	0	72	81	141	0	159	0	271	293	136	0	72	0	248	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	0	1	1	0	0	1	0	1	0	1	0	1	1	1	1	1				10	16	62.5%
2回目	測定結果	123	0	72	81	142	0	159	0	271	293	136	0	71	0	248	0	MH025463	MH025463	1				
	正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0							
	フラグ1	0	1	1	1	0	0	1	0	1	0	1	0	1	1	1	1				10	16	62.5%	
(1)ー49	1回目	測定結果	118	124	68	72	141	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	0	0	1	1	0	1	0	1	0	1	0	1	1	1	0	1				9	16	56.3%
2回目	測定結果	118	123	69	72	141	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1				
	正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247							
	フラグ1	0	0	0	1	0	1	0	1	0	1	0	1	1	1	0	1				8	16	50.0%	
(1)ー50	1回目	測定結果	135	0	69	81	128	141	159	0	271	293	132	0	72	0	244	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	0	1	0	0	1	0	0	1	0	1	0	0	1	1	1	1				8	16	50.0%
2回目	測定結果	135	0	69	81	128	141	159	0	271	292	132	0	72	0	245	0	MH025286	MH025286	1				
	正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0							
	フラグ1	0	1	0	0	1	0	0	1	0	1	0	0	1	1	0	1				7	16	43.8%	
(1)ー51	1回目	測定結果	117	123	68	81	141	0	165	0	271	0	136	0	72	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	0	0	0	1	0	0	1	0	1	0	1	0	1	1	1	1				9	16	56.3%
2回目	測定結果	118	123	69	81	142	0	165	0	271	0	136	0	72	0	248	0	MH025362	MH025362	1				
	正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0							
	フラグ1	0	0	0	0	0	1	0	1	0	1	0	1	1	1	1	1				8	16	50.0%	
(1)ー52	1回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	72	133	244	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	0	1	0	0	0	1	0	1	0	1	0	1	1	0	1	1				8	16	50.0%
2回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	72	134	244	0	ホクト白1号菌	ホクト白1号菌	1				
	正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0							
	フラグ1	0	1	0	0	0	1	0	1	0	1	0	1	0	0	1	1				7	16	43.8%	
(1)ー53	1回目	測定結果	124	0	72	81	141	0	159	0	271	293	136	0	72	0	248	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	0	1	1	1	0	0	1	0	1	0	1	0	1	1	1	1				10	16	62.5%
2回目	測定結果	124	0	73	81	141	0	159	0	271	293	136	0	72	0	248	0	MH025463	MH025463	1				
	正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0							
	フラグ1	0	1	1	0	0	0	1	0	1	0	1	0	1	1	1	1				9	16	56.3%	
(1)ー54	1回目	測定結果	123	135	69	0	128	0	152	0	267	285	145	0	92	0	248	0	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	0	0	0	0	1	1	0	1	0	1	0	1	1	1	0	1				9	16	56.3%
2回目	測定結果	123	135	69	0	128	0	152	0	267	285	145	0	92	0	248	0	MH025288	MH025288	1				
	正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
	フラグ1	0	0	0	0	1	1	1	0	1	0	1	0	1	1	0	1				9	16	56.3%	
(1)ー55	1回目	測定結果	118	135	69	81	142	0	159	0	271	285	132	136	92	0	245	0	MH025066.MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	0	0	0	0	0	1	0	1	0	1	0	1	1	0	1				6	16	37.5%
2回目	測定結果	117	135	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025066	0				
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
	フラグ1	0	0	0	0	0	1	0	1	0	1	0	0	1	1	1	1				7	16	43.8%	
(1)ー56	1回目	測定結果	117	123	68	81	141	0	165	0	271	0	136	0	72	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	0	0	0	1	0	0	1	0	1	0	1	0	1	1	1	1				9	16	56.3%
2回目	測定結果	117	123	69	81	141	0	165	0	271	0	136	0	71	0	248	0	MH025362	MH025362	1				
	正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0							
	フラグ1	0	0	0	0	0	1	0	1	0	1	0	1	1	1	1	1				8	16	50.0%	
(1)ー57	1回目	測定結果	135	0	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	1	0	0	0	1	0	1	0	1	0	0	1	1	1	1				8	16	50.0%
2回目	測定結果	135	0	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1				
	正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
	フラグ1	0	1	0	0	0	1	0	1	0	1	0	0	1	1	1	1				8	16	50.0%	
(1)ー58	1回目	測定結果	117	123	69	0	128	141	159	0	267	270	136	0	72	0	209	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	0	0	0	0	1	1	0	0	1	0	1	0	1	1	0	1				8	16	50.0%
2回目	測定結果	118	123	69	0	128	141	159	0	267	270	136	0	72	0	209	0	MH025064	MH025064	1				
	正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0							

(1)ー60	1回目	測定結果	117	0	64	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1				
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0							
		フラグ1	0	1	1	0	0	1	0	1	0	1	0	1	1	0	1	1				9		16	56.3%
	2回目	測定結果	117	0	64	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1				
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0							
		フラグ1	0	1	1	0	0	1	0	1	0	1	0	1	1	0	1	1				9		16	56.3%
		totalフラグ計	0	26	7	10	16	46	0	54	3	56	0	32	59	50	38	60		正解 品種数	100			平均 マーカ－ 一致率	48.8%
		一致率	26 11%		17 7%		62 26%		54 23%		59 25%		32 13%		109 45%		98 41%								

ブナシメジ品種識別検査結果票(B社・許容範囲±2)

サンプル 番号	検査	項目	msHm-1		msHm-2		msHm-3		msHm-4		msHm-5		msHm-6		msHm-7		msHm-8		判定品種	正解品種名	品種名一致	一致 フラグメント数	最大 フラグメント数	一致率
(1)ー1	1回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1			
		正解値 フラグ1	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247					16	16
	2回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1			
		正解値 フラグ1	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
(1)ー2	1回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1			
		正解値 フラグ1	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0					10	16
	2回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1			
		正解値 フラグ1	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0					10	16
(1)ー3	1回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1			
		正解値 フラグ1	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0					12	16
	2回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1			
		正解値 フラグ1	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
(1)ー4	1回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1			
		正解値 フラグ1	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247					12	16
	2回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1			
		正解値 フラグ1	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247					12	16
(1)ー5	1回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1			
		正解値 フラグ1	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0					12	16
	2回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1			
		正解値 フラグ1	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0					12	16
(1)ー6	1回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1			
		正解値 フラグ1	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0					12	16
	2回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1			
		正解値 フラグ1	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
(1)ー7	1回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1			
		正解値 フラグ1	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0					12	16
	2回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1			
		正解値 フラグ1	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0					12	16
(1)ー8	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0			
		正解値 フラグ1	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0					12	16
	2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0			
		正解値 フラグ1	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0					12	16
(1)ー9	1回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1			
		正解値 フラグ1	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0					12	16
	2回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1			
		正解値 フラグ1	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
(1)ー10	1回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1			
		正解値 フラグ1	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0					10	16
	2回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1			
		正解値 フラグ1	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0					10	16
(1)ー11	1回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1			
		正解値 フラグ1	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0					12	16
	2回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1			
		正解値 フラグ1	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0					12	16

(1)ー12	1回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー13	2回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー14	1回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー15	2回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー16	1回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー17	2回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー18	1回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー19	2回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー20	1回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー21	2回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー22	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー23	2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー24	1回目	測定結果	117	133	66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				10	16	62.5%
(1)ー25	2回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				10	16	62.5%
(1)ー26	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー27	2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー28	1回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー29	2回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー30	1回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー31	2回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー32	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー33	2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%

(1)ー24	1回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1				12	16	75.0%
(1)ー25	2回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1				12	16	75.0%
(1)ー26	1回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー27	2回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー28	1回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー29	2回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー30	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー31	2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー32	1回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー33	2回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー34	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー35	2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー36	1回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー37	2回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー38	1回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー39	2回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー40	1回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー41	2回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%

(1)ー36	1回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1					
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0					12	16	75.0%	
		フラグ1		1		1	1	1		1		1	1	1	1		1	1								
2回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1						
	正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0									
	フラグ1		1		1	1	1		1		1	1	1	1		1	1					12	16	75.0%		
(1)ー37	1回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1					
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0					12	16	75.0%	
		フラグ1		1		1	1	1		1		1	1	1		0	1	1								
2回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1						
	正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						12	16	75.0%	
	フラグ1		1		1	1	1		1		1	1	1		0	1	1					12	16	75.0%		
(1)ー38	1回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1					
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						12	16	75.0%
		フラグ1		1		1	1	1		1		1	1	1	1		0	1					12	16	75.0%	
2回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1						
	正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							12	16	75.0%
	フラグ1		1		1	1	1		1		1	1	1		0	1	1						12	16	75.0%	
(1)ー39	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0					
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						12	16	75.0%
		フラグ1		1		1	1	1		1		1	1	1	1		1	1					12	16	75.0%	
2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0						
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							12	16	75.0%
	フラグ1		1		1	1	1		1		1	1	1	1		1	1						12	16	75.0%	
(1)ー40	1回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1					
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						12	16	75.0%
		フラグ1		1		1	1	1		1		1	1	1	1		1	1						12	16	75.0%
2回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1						
	正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							12	16	75.0%
	フラグ1		1		1	1	1		1		1	1	1	1		1	1						12	16	75.0%	
(1)ー41	1回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1					
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						12	16	75.0%
		フラグ1		1		1	1	1		1		1	1	1	1		0	1					12	16	75.0%	
2回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1						
	正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							12	16	75.0%
	フラグ1		1		1	1	1		1		1	1	1	1		0	1						12	16	75.0%	
(1)ー42	1回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1					
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						12	16	75.0%
		フラグ1		1		1	1	1	1		1	1	1	1	1		1	1					12	16	75.0%	
2回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1						
	正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						12	16	75.0%	
	フラグ1		1		1	1	1	1		1	1	1	1	1		1	1					12	16	75.0%		
(1)ー43	1回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1					
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						12	16	75.0%
		フラグ1		1		1	1	1		1		1	1	1	1		1	1					12	16	75.0%	
2回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1						
	正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247							12	16	75.0%
	フラグ1		1		1	1	1		1		1	1	1	1		1	1						12	16	75.0%	
(1)ー44	1回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1					
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						12	16	75.0%
		フラグ1		1		1	1	1		1		1	1	1	1		1	1					12	16	75.0%	
2回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1						
	正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0							12	16	75.0%
	フラグ1		1		1	1	1		1		1	1	1	1		1	1						12	16	75.0%	
(1)ー45	1回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1					
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						12	16	75.0%
		フラグ1		1		1	1	1	1		1	1	1	1	1		1	1					12	16	75.0%	
2回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1						
	正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0							12	16	75.0%
	フラグ1		1		1	1	1	1		1	1	1	1	1		1	1						12	16	75.0%	
(1)ー46	1回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1					
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						10	16	62.5%
		フラグ1		1		1	1	1	1	0		1	0	1	1		1	1					10	16	62.5%	
2回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1						
	正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0									

(1)ー48	1回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1				
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0					12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1							
(1)ー49	2回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1				
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0					12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1							
(1)ー50	1回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1				
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247					12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1							
(1)ー51	2回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1				
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247					12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1							
(1)ー52	1回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1				
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0					12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1							
(1)ー53	2回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1				
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0					12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1							
(1)ー54	1回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1				
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0					12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1							
(1)ー55	2回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1				
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0					12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1							
(1)ー56	1回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1				
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0					12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1							
(1)ー57	2回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1				
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0					12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1							
(1)ー58	1回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1				
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0					12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1							
(1)ー59	2回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1				
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0					12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1							
(1)ー60	1回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1				
		正解値	133	0	67	79	140	0	155	0	269	284	130	134	92	0	243	0					6	16	37.5%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1							
(1)ー61	2回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1				
		正解値	133	0	67	79	140	0	155	0	269	284	130	134	92	0	243	0					6	16	37.5%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1							
(1)ー62	1回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1				
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0					12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1							
(1)ー63	2回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1				
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0					12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1							
(1)ー64	1回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1				
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0					10	16	62.5%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1							
(1)ー65	2回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1				
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0					10	16	62.5%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1							

(1)ー60	1回目	測定結果	117	62	79	145	150	270	284	131	70	243	MH025397	MH025397	1				
		正解値 フラグ1	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0	
	2回目	測定結果	117	62	79	145	150	270	284	131	70	243	MH025397	MH025397	1		12	16	75.0%
		正解値 フラグ1	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0	
	totalフラグ計		60	60	60	60	60	60	60	54	60	60	54	60	60	50	56	60	
	一致率		50%	50%	50%	48%	50%	48%	46%	48%									
														正解 品種数	100		平均 マーカ ー一致率	73.5%	

ブナシメジ品種識別検査結果票(F社・許容範囲±2)

サンプル 番号	検査	項目	msHm-1		msHm-2		msHm-3		msHm-4		msHm-5		msHm-6		msHm-7		msHm-8		判定品種	正解品種名	品種名一致	一致 フラグメント数	最大 フラグメント数	一致率
(1)ー1	1回目	測定結果	117	122	66	71	139	0	161	0	270	281	131	135	70	0	217	247	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1	1				15	16	93.8%
(1)ー2	2回目	測定結果	117	122	66	71	139	0	161	0	270	281	135	0	70	0	217	246	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
(1)ー3	1回目	測定結果	122	146	66	71	127	139	149	162	270	0	131	135	70	0	217	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	0	1	1	1	1	1	0	1	1	1	1	1	1	1	1				14	16	87.5%
(1)ー4	2回目	測定結果	122	146	66	71	127	139	149	161	270	0	131	135	70	0	217	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	0	1	1	1	1	1	0	1	1	1	1	1	1	1	1				14	16	87.5%
(1)ー5	1回目	測定結果	122	134	66	0	127	145	149	155	266	284	131	0	70	90	217	247	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	1	1	1	1	1	0	1	0	1	1	0	1	0	0	0	0				10	16	62.5%
(1)ー6	2回目	測定結果	122	134	66	0	127	0	149	0	266	284	131	143	70	90	×	0	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	0	1	0	0	1	1				13	16	81.3%
(1)ー7	1回目	測定結果	117	122	66	70	139	0	162	0	270	281	131	135	70	0	217	247	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1	1				15	16	93.8%
(1)ー8	2回目	測定結果	117	122	66	71	139	0	161	0	270	281	135	0	70	0	217	247	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
(1)ー9	1回目	測定結果	117	122	66	79	127	139	162	0	270	0	135	0	70	0	247	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー10	2回目	測定結果	117	122	66	79	139	0	161	0	270	0	135	0	70	0	246	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
(1)ー11	1回目	測定結果	134	0	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
(1)ー12	2回目	測定結果	134	0	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
(1)ー13	1回目	測定結果	117	0	×	0	146	0	149	0	270	284	131	0	×	0	243	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	1	1	0	0	0	1	1	1	1	1	1	1	0	0	1	1				12	16	75.0%
(1)ー14	2回目	測定結果	117	0	62	79	146	0	149	0	270	284	131	0	70	0	243	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1				15	16	93.8%
(1)ー15	1回目	測定結果	117	134	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025066・MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
(1)ー16	2回目	測定結果	117	134	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025066・MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
(1)ー17	1回目	測定結果	134	146	66	79	127	139	156	0	270	292	131	0	70	0	243	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	0	1	1	1	1				14	16	87.5%
(1)ー18	2回目	測定結果	134	146	67	79	127	139	155	0	270	292	131	0	70	0	242	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	0	1	1	1	1	1	1	1	1	0	1	1	1	1	1				14	16	87.5%
(1)ー19	1回目	測定結果	122	146	66	71	127	139	149	161	270	0	131	135	70	0	217	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	0	1	1	1	1	1	0	1	1	1	1	1	1	1	1				14	16	87.5%
(1)ー20	2回目	測定結果	117	122	66	71	127	139	149	161	270	0	131	135	70	0	217	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	0	1	1	1	1	1	0	1	1	1	1	1	1	1	1				14	16	87.5%
(1)ー21	1回目	測定結果	117	122	66	0	127	139	156	0	266	270	131	135	70	0	206	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	0	1	1	1	1	1				15	16	93.8%
(1)ー22	2回目	測定結果	117	122	66	0	127	139	156	0	266	270	135	143	70	0	×	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	0	1	1	0	1	1				14	16	87.5%

(1)ー12	1回目	測定結果	134	0	66	79	145	0	156	0	270	0	131	0	70	132	243	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1						
	2回目	測定結果	134	0	66	79	145	0	155	0	270	0	131	0	×	0	243	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1						
(1)ー13	1回目	測定結果	116	122	66	79	139	0	161	0	270	0	135	0	70	0	246	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0				16	16	100.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1						
	2回目	測定結果	117	122	66	79	139	0	161	0	270	0	135	0	70	0	246	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0				16	16	100.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1						
(1)ー14	1回目	測定結果	122	146	71	79	139	0	156	0	270	292	131	135	70	0	247	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0				14	16	87.5%
		フラグ1	1	0	1	1	1	1	1	1	1	1	0	1	1	1	1	1						
	2回目	測定結果	122	146	71	79	139	0	155	0	×	0	135	0	70	0	247	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	1	0	1	1	1	1	1	1	1	0	0	1	1	1	1	1				13	16	81.3%
(1)ー15	1回目	測定結果	134	0	66	79	145	0	155	0	270	0	131	0	70	132	243	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1						
	2回目	測定結果	134	0	66	79	145	0	155	0	270	0	131	0	70	132	242	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1						
(1)ー16	1回目	測定結果	134	0	66	79	127	139	155	0	270	292	131	0	70	0	243	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	1	1	1	0	1	1	1	1	1						
	2回目	測定結果	134	146	67	79	127	139	156	0	270	292	131	135	70	0	243	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0				15	16	93.8%
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	1	1	1	1	1						
(1)ー17	1回目	測定結果	×	0	66	79	139	0	155	0	270	284	131	135	90	0	243	0	MH025066・MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				14	16	87.5%
		フラグ1	0	0	1	1	1	1	1	1	1	1	1	1	1	1	1	1						
	2回目	測定結果	117	134	66	79	139	0	155	0	270	284	131	135	90	0	243	0	MH025066・MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				16	16	100.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1						
(1)ー18	1回目	測定結果	122	146	66	71	127	139	149	161	270	0	131	135	70	0	217	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0				14	16	87.5%
		フラグ1	1	0	1	1	1	1	1	0	1	1	1	1	1	1	1	1						
	2回目	測定結果	122	146	66	71	127	139	149	161	270	0	131	135	70	0	217	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0				14	16	87.5%
		フラグ1	1	0	1	1	1	1	1	0	1	1	1	1	1	1	1	1						
(1)ー19	1回目	測定結果	117	134	66	79	139	0	155	0	270	284	131	135	90	0	242	0	MH025066・MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				16	16	100.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1						
	2回目	測定結果	117	134	66	79	139	0	156	0	270	284	131	135	90	0	242	0	MH025066・MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				16	16	100.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1						
(1)ー20	1回目	測定結果	117	134	66	79	139	0	155	0	270	284	131	135	70	90	246	0	MH025066・MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				10	16	62.5%
		フラグ1	1	0	1	1	1	1	1	0	1	1	0	1	0	1	0	1						
	2回目	測定結果	117	134	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025066・MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				10	16	62.5%
		フラグ1	1	0	1	1	1	1	1	0	1	1	0	1	0	1	0	1						
(1)ー21	1回目	測定結果	117	122	66	71	139	145	156	161	270	281	131	0	70	0	217	246	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				8	16	50.0%
		フラグ1	1	0	1	0	1	0	0	1	1	0	1	0	1	1	0	0						
	2回目	測定結果	117	122	66	71	139	0	161	0	270	281	135	0	70	132	217	247	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				8	16	50.0%
		フラグ1	1	0	1	0	1	1	0	1	1	0	1	0	0	1	0	0						
(1)ー22	1回目	測定結果	134	0	66	79	145	0	156	0	270	0	131	0	70	132	243	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0				12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	0	1	0	0	1	1	1						
	2回目	測定結果	134	0	66	79	145	0	155	0	270	0	131	0	70	132	243	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269													

(1)ー24	1回目	測定結果	×	0	66	0	127	0	149	156	266	284	131	0	70	90	221	0	MH025288	MH025288	1				
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
		フラグ1	0	0	1	1	1	1	1	1	0	1	0	1	1	0	1	1					11	16	68.8%
		測定結果	122	134	66	0	127	0	149	156	266	284	131	143	70	90	×	0	MH025288	MH025288	1				
(1)ー25	2回目	正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
		フラグ1	1	1	1	1	1	1	1	0	1	1	1	0	1	0	0	1					12	16	75.0%
		測定結果	117	122	66	0	127	139	155	0	266	270	131	135	70	0	246	0	MH025064	MH025064	1				
		正解値	115	121	67	0	127	139	155	0	265	269	134	72	72	0	207	0							
(1)ー26	1回目	フラグ1	1	0	0	1	0	0	1	1	0	1	0	0	1	0	1	0					9	16	56.3%
		測定結果	117	122	66	0	127	139	155	0	×	0	135	0	70	0	206	0	MH025064	MH025064	1				
		正解値	115	121	67	0	127	139	155	0	265	269	134	72	72	0	207	0							
		フラグ1	1	0	1	0	0	1	1	1	0	1	1	0	0	1	0	1					9	16	56.3%
(1)ー27	2回目	測定結果	134	0	66	79	139	0	155	0	270	284	131	135	90	0	243	0	MH025528	MH025528	1				
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0							
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1					15	16	93.8%
		測定結果	134	0	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025528	MH025528	1				
(1)ー28	1回目	正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0							
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1					15	16	93.8%
		測定結果	117	134	66	79	139	145	155	0	270	284	131	135	70	0	242	0	MH025066・MH025012	MH025066	0				
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
(1)ー29	2回目	フラグ1	1	0	0	1	0	1	0	1	0	1	0	0	1	0	0	0					7	16	43.8%
		測定結果	117	134	66	79	139	0	156	0	270	284	131	135	90	0	242	0	MH025066・MH025012	MH025066	0				
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	1	0	1	0	1	1	0	1	0	1	0	0	0	1	0	1					7	16	43.8%
(1)ー30	1回目	測定結果	116	122	66	79	139	0	161	0	270	0	135	0	70	0	246	0	MH025362	MH025362	1				
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0							
		フラグ1	1	0	1	1	1	1	1	0	1	1	0	1	0	1	0	1					11	16	68.8%
		測定結果	116	122	66	79	139	0	161	0	270	0	135	0	70	0	247	0	MH025362	MH025362	1				
(1)ー31	2回目	正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0							
		フラグ1	1	0	1	1	1	1	1	0	1	0	1	0	0	1	0	1					10	16	62.5%
		測定結果	117	134	66	79	139	145	155	0	270	284	131	135	70	90	243	0	MH025066・MH025012	MH025012	0				
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
(1)ー32	1回目	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1					14	16	87.5%
		測定結果	117	134	66	79	139	0	155	0	270	284	131	135	90	0	242	0	MH025066・MH025012	MH025012	0				
		正解値	115	133	67	79																			

(1)－36	1回目	測定結果	134	0	66	79	145	0	156	0	270	284	131	135	90	0	242	0	MH025528	MH025528	1					
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0					15	16	93.8%	
		フラグ1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1	1								
2回目	測定結果	134	0	66	79	139	0	156	0	270	284	131	135	90	0	242	0	MH025528	MH025528	1						
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0								
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1					16	16	100.0%	
(1)－37	1回目	測定結果	117	0	62	79	146	0	149	0	270	284	131	0	70	0	242	0	MH025397	MH025397	1					
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0								
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1					15	16	93.8%	
2回目	測定結果	116	0	62	79	146	0	149	0	270	284	131	0	70	0	243	0	MH025397	MH025397	1						
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0								
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1					15	16	93.8%	
(1)－38	1回目	測定結果	117	122	66	0	127	0	155	0	266	284	131	0	90	0	×	0	MH025288	MH025288	1					
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0								
		フラグ1	1	0	1	1	1	1	1	0	1	1	0	1	1	1	0	1					12	16	75.0%	
2回目	測定結果	122	134	66	0	127	0	149	0	266	284	131	143	70	90	243	0	MH025288	MH025288	1						
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0								
		フラグ1	1	1	1	1	1	1	1	1	1	1	0	1	0	0	1									
(1)－39	1回目	測定結果	111	117	66	79	139	145	156	0	270	284	131	135	90	0	242	0	MH025066・MH025012	MH025012	0					
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0								
		フラグ1	1	0	1	1	1	0	1	1	1	1	1	1	1	1	1	1					14	16	87.5%	
2回目	測定結果	117	134	66	79	139	0	155	0	270	284	131	135	90	0	243	0	MH025066・MH025012	MH025012	0						
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0								
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1					16	16	100.0%	
(1)－40	1回目	測定結果	134	0	66	79	139	0	156	0	270	284	135	0	70	90	242	0	MH025528	MH025528	1					
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0								
		フラグ1	1	1	1	1	1	1	1	1	1	1	0	1	1	0	1	1					14	16	87.5%	
2回目	測定結果	134	0	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025528	MH025528	1						
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0								
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1					16	16	100.0%	
(1)－41	1回目	測定結果	122	134	66	0	128	0	149	155	266	284	131	0	70	90	206	221	MH025288	MH025288	1					
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0								
		フラグ1	1	1	1	1	1	1	1	0	1	1	0	1	1	0	1	0					12	16	75.0%	
2回目	測定結果	122	134	66	0	127	0	149	155	266	284	143	0	70	90	×	0	MH025288	MH025288	1						
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0								
		フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	1	0	0	1					13	16	81.3%	
(1)－42	1回目	測定結果	134	0	66	79	127	139	156	0	270	292	131	0	70	0	242									

(1)ー48	1回目	測定結果	122	146	71	79	139	0	156	0	270	292	135	0	70	0	246	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0				15	16	93.8%
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	1	1	1	1	1						
	2回目	測定結果	122	146	71	79	139	0	155	0	270	292	135	0	70	0	246	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0				15	16	93.8%
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	1	1	1	1	1						
(1)ー49	1回目	測定結果	117	122	66	71	139	0	161	0	270	281	135	0	70	0	217	246	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				16	16	100.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1						
	2回目	測定結果	117	122	66	71	139	0	161	0	270	281	135	0	70	0	217	246	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				16	16	100.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1						
(1)ー50	1回目	測定結果	134	0	66	79	127	139	155	0	270	292	131	135	70	0	242	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0				16	16	100.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1						
	2回目	測定結果	134	146	67	79	127	139	156	0	270	292	131	0	70	0	217	242	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	0	1	1	1	0				13	16	81.3%
(1)ー51	1回目	測定結果	117	122	66	79	139	0	161	0	270	0	135	0	70	0	246	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0				16	16	100.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1						
	2回目	測定結果	116	122	66	79	139	0	161	0	270	0	135	0	70	0	246	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0				16	16	100.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1						
(1)ー52	1回目	測定結果	122	134	66	79	145	0	156	0	270	0	131	135	70	132	242	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0				13	16	81.3%
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	0	1	0	1	1						
	2回目	測定結果	122	134	66	79	145	0	156	0	270	0	131	135	70	132	242	246	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0				12	16	75.0%
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	0	1	0	1	0						
(1)ー53	1回目	測定結果	122	146	71	79	139	0	156	0	270	292	135	0	70	0	246	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0				15	16	93.8%
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	1	1	1	1	1						
	2回目	測定結果	122	146	71	79	127	139	155	0	270	292	135	0	70	132	246	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0				13	16	81.3%
		フラグ1	1	0	1	1	1	0	1	1	1	1	1	1	1	0	1	1						
(1)ー54	1回目	測定結果	122	134	66	0	127	0	149	168	266	284	131	143	70	90	217	246	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0				11	16	68.8%
		フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	0	1	0	0						
	2回目	測定結果	122	134	66	0	127	145	149	155	266	284	131	143	70	90	217	0	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0				11	16	68.8%
		フラグ1	1	1	1	1	1	0	1	0	1	1	1	0	1	0	0	1						
(1)ー55	1回目	測定結果	117	134	66	79	139	0	156	0	270	284	131	135	90	0	217	242	MH025066・MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0						
	2回目	測定結果	117	134	66	79	139	0	155	0	270	284	131	135	90	0	242	0	MH025066・MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				16	16	100.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1						
(1)ー56	1回目	測定結果	117	122	66	79	139	0	161	0	270	0	131	135	70	0	246	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1	1						
	2回目	測定結果	117	122	66	79	139	0	161	0	270	0	135	0	70	0	246	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0				16	16	100.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1						
(1)ー57	1回目	測定結果	134	0	66	79	139	0	155	0	270	284	131	135	90	0	217	242	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0						
	2回目	測定結果	134	0	66	79	139	0	155	0	270	284	131	135	90	0	217	242	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0						
(1)ー58	1回目	測定結果	117	122	66	0	127	139	156	0	266	270	135	0	70	0	217	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1						
	2回目	測定結果	117	122	67	0	127	139	156	0	266	270	135	0	70	0	217	0	MH025064</					

(1)ー60	1回目	測定結果	117	0	62	79	127	146	149	0	270	284	131	0	70	0	242	0	MH025397	MH025397	1					
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0								
		フラグ1	1	1	1	1	1	0	1	1	1	1	1	1	1	0	1	1					14		16	87.5%
	2回目	測定結果	117	0	62	79	146	0	149	0	270	284	131	0	70	0	242	0	MH025397	MH025397	1					
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0								
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1					15		16	93.8%
	totalフラグ計		58	35	59	51	56	52	53	48	56	47	54	41	46	44	45	55								
			93		110		108		101		103		95		90		100									
		一致率	39%		46%		45%		42%		43%		40%		38%		42%									
																				正解 品種数	100				平均 マーカ ー致率	87.1%

ブナシメジ品種識別検査結果票(L社・許容範囲±2)

サンプル 番号	検査	項目	msHm-1		msHm-2		msHm-3		msHm-4		msHm-5		msHm-6		msHm-7		msHm-8		判定品種	正解品種名	品種名一致	一致 フラグメント数	最大 フラグメント数	一致率
(1)-1	1回目	測定結果	117	123	69	73	141	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)-2	2回目	測定結果	118	124	69	73	142	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	0	0	1	1	0	1	0	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)-2	1回目	測定結果	123	0	69	73	128	142	152	165	270	0	136	0	72	0	219	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	1	1	1	1	0	0	0	1	1	0	1	1	1	1	1				12	16	75.0%
(1)-3	2回目	測定結果	124	0	69	73	128	141	152	165	271	0	136	0	72	0	219	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	0	1	1	1	1	1	0	0	1	1	0	1	1	1	1	1				12	16	75.0%
(1)-3	1回目	測定結果	123	135	69	0	128	0	152	0	267	285	145	0	92	0	0	0	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	1	1	1	1	1	1	0	1	1	1	0	1	1	1	0	1				13	16	81.3%
(1)-4	2回目	測定結果	123	135	69	0	128	0	152	0	267	285	146	0	92	0	0	0	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	1	1	1	1	1	1	0	1	1	1	0	1	1	1	0	1				13	16	81.3%
(1)-4	1回目	測定結果	118	124	69	73	141	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	0	0	1	1	1	1	0	1	1	1	1	1	1	1	1	1				13	16	81.3%
(1)-5	2回目	測定結果	118	124	69	73	141	0	165	0	271	282	136	0	71	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	0	0	1	1	1	1	0	1	1	1	1	1	1	1	1	1				13	16	81.3%
(1)-5	1回目	測定結果	117	124	69	81	141	0	165	0	271	0	137	0	72	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	0	1	1	1	0	1	1	0	1	1	1	1	1	1				13	16	81.3%
(1)-6	2回目	測定結果	118	124	69	81	141	0	165	0	271	0	136	0	71	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	0	0	1	1	1	1	0	1	1	1	1	1	1	1	1	1				13	16	81.3%
(1)-6	1回目	測定結果	135	0	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)-7	2回目	測定結果	135	0	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1	1				14	16	87.5%
(1)-7	1回目	測定結果	118	0	65	81	149	0	152	0	271	285	132	0	72	0	245	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	0	1	1	1	0	1	0	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)-8	2回目	測定結果	118	0	65	81	149	0	153	0	271	285	132	0	71	0	244	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	0	1	1	1	0	1	0	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)-8	1回目	測定結果	118	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066,MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				14	16	87.5%
(1)-9	2回目	測定結果	118	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066,MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				14	16	87.5%
(1)-9	1回目	測定結果	135	0	69	81	128	141	159	0	271	293	132	0	72	0	244	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)-10	2回目	測定結果	135	0	69	81	129	141	159	0	271	293	132	0	72	0	244	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)-10	1回目	測定結果	123	0	69	73	128	141	152	165	271	0	136	0	72	0	219	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	1	1	1	1	1	0	0	1	1	0	1	1	1	1	1				13	16	81.3%
(1)-11	2回目	測定結果	123	0	69	73	128	142	152	165	270	0	136	0	71	0	219	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	1	1	1	1	1	0	0	1	1	0	1	1	1	1	1				12	16	75.0%
(1)-11	1回目	測定結果	117	123	69	0	128	141	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)-11	2回目	測定結果	118	123	69	0	128	142	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	0	1	1	1	1	0	0	1	1	1	1	1	1	1	1	1				13	16	81.3%

(1)ー12	1回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	72	133	244	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	0	1	0	1	1	1	1	1	0	1	1				13	16	81.3%
2回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	72	133	244	0	ホクト白1号菌	ホクト白1号菌	1				
	正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0							
	フラグ1	1	1	1	1	1	0	1	0	1	1	1	1	1	0	1	1				13	16	81.3%	
(1)ー13	1回目	測定結果	117	124	69	81	141	0	165	0	271	0	136	0	72	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	0	1	1	1	1	1	0	1	1	1	1	1	1	1	1				14	16	87.5%
2回目	測定結果	118	123	69	81	141	0	165	0	271	0	136	0	71	0	248	0	MH025362	MH025362	1				
	正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0							
	フラグ1	0	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1				14	16	87.5%	
(1)ー14	1回目	測定結果	123	0	73	81	141	0	159	0	271	293	136	0	72	0	248	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1				15	16	93.8%
2回目	測定結果	124	0	73	81	141	0	159	0	271	293	136	0	71	0	248	0	MH025463	MH025463	1				
	正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0							
	フラグ1	0	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1				14	16	87.5%	
(1)ー15	1回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	72	133	244	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	0	1	0	1	1	1	1	1	0	1	1				13	16	81.3%
2回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	72	133	244	0	ホクト白1号菌	ホクト白1号菌	1				
	正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0							
	フラグ1	1	1	1	1	1	0	1	0	1	1	1	1	1	0	1	1				13	16	81.3%	
(1)ー16	1回目	測定結果	135	0	69	81	128	141	159	0	271	293	132	0	72	0	244	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1				15	16	93.8%
2回目	測定結果	135	0	69	81	129	141	159	0	271	293	132	0	72	0	244	0	MH025286	MH025286	1				
	正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0							
	フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1				15	16	93.8%	
(1)ー17	1回目	測定結果	117	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066,MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1				15	16	93.8%
2回目	測定結果	118	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066,MH025012	MH025066	0				
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
	フラグ1	0	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1				14	16	87.5%	
(1)ー18	1回目	測定結果	123	0	69	72	128	141	152	165	271	0	136	0	72	0	219	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	1	1	1	1	1	1	0	0	1	1	0	1	1	1	1				13	16	81.3%
2回目	測定結果	123	0	69	72	128	141	152	165	271	0	136	0	71	0	219	0	MH025462	MH025462	1				
	正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0							
	フラグ1	1	1	1	1	1	1	1	0	0	1	1	0	1	1	1	1				13	16	81.3%	
(1)ー19	1回目	測定結果	117	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066,MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1				15	16	93.8%
2回目	測定結果	118	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066,MH025012	MH025012	0				
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
	フラグ1	0	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1				14	16	87.5%	
(1)ー20	1回目	測定結果	117	135	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025066,MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	0	1	0	1	1	1	1	1	1	1				14	16	87.5%
2回目	測定結果	118	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066,MH025012	MH025012	0				
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
	フラグ1	0	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1				14	16	87.5%	
(1)ー21	1回目	測定結果	118	123	69	73	142	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	0	1	1	1	0	1	0	1	1	1	1	1	1	1	1	1				13	16	81.3%
2回目	測定結果	118	123	69	73	141	0	165	0	271	282	136	0	71	0	219	248	MH025369	MH025369	1				
	正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247							
	フラグ1	0	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1				14	16	87.5%	
(1)ー22	1回目	測定結果	135	0	68	81	148	0	159	0	271	0	132	0	72	133	244	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	0	1	0	1	1	1	1	1	0	1	1				13	16	81.3%
2回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	71	134	244	0	ホクト白1号菌	ホクト白1号菌	1				
	正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0							
	フラグ1	1	1	1	1	1	0	1	0	1														

(1)ー24	1回目	測定結果	123	135	69	0	128	0	152	0	267	285	145	0	92	0	244	0	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	1	1	1	1	1	1	1	0	1	1	1	0	1	1	1	0	1				13	16
2回目	測定結果	124	135	68	0	128	0	152	0	267	285	145	0	92	0	244	0	MH025288	MH025288	1				
	正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
	フラグ1	0	1	1	1	1	1	1	0	1	1	1	0	1	1	1	0	1				12	16	75.0%
(1)ー25	1回目	測定結果	117	123	69	0	128	141	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				15	16
2回目	測定結果	117	123	68	0	128	141	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1				
	正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0							
	フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー26	1回目	測定結果	135	0	68	81	141	0	159	0	271	285	132	136	92	0	245	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				15	16
2回目	測定結果	135	0	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1				
	正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
	フラグ1	1	1	1	1	0	1	0	1	1	1	1	1	1	1	1	1	1				14	16	87.5%
(1)ー27	1回目	測定結果	117	135	68	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025066,MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	0	1	0	1	1	1	1	1	1	1	1	1	1				14	16
2回目	測定結果	117	135	69	81	142	0	159	0	270	285	132	136	92	0	244	0	MH025066,MH025012	MH025066	0				
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
	フラグ1	1	1	1	1	0	1	0	1	1	1	1	1	1	1	1	1	1				14	16	87.5%
(1)ー28	1回目	測定結果	117	123	68	81	142	0	165	0	271	0	136	0	72	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	0	1	0	1	1	1	1	1	1	1	1	1	1				14	16
2回目	測定結果	118	123	69	81	141	0	165	0	271	0	136	0	71	0	248	0	MH025362	MH025362	1				
	正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0							
	フラグ1	0	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				14	16	87.5%
(1)ー29	1回目	測定結果	117	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066,MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				15	16
2回目	測定結果	118	135	69	81	142	0	159	0	271	285	132	137	92	0	244	0	MH025066,MH025012	MH025012	0				
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
	フラグ1	0	1	1	1	0	1	0	1	1	1	1	1	1	1	1	1	1				13	16	81.3%
(1)ー30	1回目	測定結果	118	0	64	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	0	1	1	1	0	1	0	1	1	1	1	1	0	1	1	1	1				12	16
2回目	測定結果	118	0	65	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1				
	正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0							
	フラグ1	0	1	1	1	0	1	0	1	1	1	1	1	0	1	1	1	1				12	16	75.0%
(1)ー31	1回目	測定結果	117	123	69	0	128	141	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				15	16
2回目	測定結果	117	123	69	0	128	141	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1				
	正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0							
	フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー32	1回目	測定結果	123	0	72	81	141	0	159	0	271	293	136	0	72	0	248	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				15	16
2回目	測定結果	123	0	73	81	141	0	159	0	271	293	136	0	71	0	248	0	MH025463	MH025463	1				
	正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0							
	フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー33	1回目	測定結果	117	0	65	81	149	0	152	0	271	285	132	0	72	0	245	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	1	1	1	1	0	1	0	1	1	1	1	1	1	0	1	1	1				13	16
2回目	測定結果	117	0	65	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1				
	正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0							
	フラグ1	1	1	1	1	0	1	0	1	1	1	1	1	1	0	1	1	1				13	16	81.3%
(1)ー34	1回目	測定結果	117	123	69	0	128	141	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				15	16
2回目	測定結果	117	123	69	0	128	141	159	0	267	271	136	0	71	0	209	0	MH025064	MH025064	1				
	正解値	115																						

(1)ー36	1回目	測定結果	135	0	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1				
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0					15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1							
2回目	測定結果	135	0	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1					
	正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0								
	フラグ1	1	1	1	1	0	1	0	1	1	1	1	1	1	1	1	1					14	16	87.5%	
(1)ー37	1回目	測定結果	117	0	64	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1				
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0							
		フラグ1	1	1	1	1	0	1	0	1	1	1	1	1	0	1	1	1					13	16	81.3%
2回目	測定結果	117	0	65	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1					
	正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0								
	フラグ1	1	1	1	1	0	1	0	1	1	1	1	1	0	1	1	1					13	16	81.3%	
(1)ー38	1回目	測定結果	123	135	69	0	128	0	152	0	267	285	145	0	92	0	245	0	MH025288	MH025288	1				
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
		フラグ1	1	1	1	1	1	1	0	1	1	1	0	1	1	1	0	1					13	16	81.3%
2回目	測定結果	123	135	69	0	128	0	152	0	267	285	145	0	92	0	244	0	MH025288	MH025288	1					
	正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0								
	フラグ1	1	1	1	1	1	1	0	1	1	1	0	1	1	1	0	1					13	16	81.3%	
(1)ー39	1回目	測定結果	117	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025012	0				
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1					15	16	93.8%
2回目	測定結果	117	135	69	81	141	0	159	0	271	285	132	136	92	0	245	0	MH025066.MH025012	MH025012	0					
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0								
	フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1					15	16	93.8%	
(1)ー40	1回目	測定結果	135	0	69	81	141	0	159	0	271	285	132	136	92	0	245	0	MH025528	MH025528	1				
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1					15	16	93.8%
2回目	測定結果	135	0	69	81	141	0	159	0	271	285	132	136	92	0	245	0	MH025528	MH025528	1					
	正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0								
	フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1					15	16	93.8%	
(1)ー41	1回目	測定結果	123	135	69	0	129	0	152	0	267	285	145	0	92	0	244	0	MH025288	MH025288	1				
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
		フラグ1	1	1	1	1	1	1	0	1	1	1	0	1	1	1	0	1					13	16	81.3%
2回目	測定結果	123	135	69	0	129	0	152	0	267	285	145	0	92	0	244	0	MH025288	MH025288	1					
	正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0								
	フラグ1	1	1	1	1	1	1	0	1	1	1	0	1	1	1	0	1					13	16	81.3%	
(1)ー42	1回目	測定結果	135	0	69	81	128	141	159	0	271	293	132	0	72	0	244	0	MH025286	MH025286	1				
		正解値	133	0	67	79	1																		

(1)ー48	1回目	測定結果	123	0	72	81	141	0	159	0	271	293	136	0	72	0	248	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1						
2回目	測定結果	123	0	72	81	142	0	159	0	271	293	136	0	71	0	248	0	MH025463	MH025463	1				
	正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0							
	フラグ1	1	1	1	1	1	0	1	0	1	1	1	1	1	1	1	1				14	16	87.5%	
(1)ー49	1回目	測定結果	118	124	68	72	141	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				13	16	81.3%
		フラグ1	0	0	1	1	1	1	0	1	1	1	1	1	1	1	1	1						
2回目	測定結果	118	123	69	72	141	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1				
	正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				14	16	87.5%	
	フラグ1	0	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1							
(1)ー50	1回目	測定結果	135	0	69	81	128	141	159	0	271	293	132	0	72	0	244	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1						
2回目	測定結果	135	0	69	81	128	141	159	0	271	292	132	0	72	0	245	0	MH025286	MH025286	1				
	正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0							
	フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				15	16	93.8%	
(1)ー51	1回目	測定結果	117	123	68	81	141	0	165	0	271	0	136	0	72	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1						
2回目	測定結果	118	123	69	81	142	0	165	0	271	0	136	0	72	0	248	0	MH025362	MH025362	1				
	正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0							
	フラグ1	0	1	1	1	1	0	1	0	1	1	1	1	1	1	1	1				13	16	81.3%	
(1)ー52	1回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	72	133	244	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	0	1	0	1	1	1	1	1	0	1	1				13	16	81.3%
2回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	72	134	244	0	ホクト白1号菌	ホクト白1号菌	1				
	正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0							
	フラグ1	1	1	1	1	1	0	1	0	1	1	1	1	1	0	1	1				13	16	81.3%	
(1)ー53	1回目	測定結果	124	0	72	81	141	0	159	0	271	293	136	0	72	0	248	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0				14	16	87.5%
		フラグ1	0	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1						
2回目	測定結果	124	0	73	81	141	0	159	0	271	293	136	0	72	0	248	0	MH025463	MH025463	1				
	正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0							
	フラグ1	0	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				14	16	87.5%	
(1)ー54	1回目	測定結果	123	135	69	0	128	0	152	0	267	285	145	0	92	0	248	0	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0				13	16	81.3%
		フラグ1	1	1	1	1	1	1	0	1	1	1	0	1	1	1	0	1						
2回目	測定結果	123	135	69	0	128	0	152	0	267	285	145	0	92	0	248	0	MH025288	MH025288	1				
	正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
	フラグ1	1	1	1	1	1	1	0	1	1	1	0	1	1	1	0	1				13	16	81.3%	
(1)ー55	1回目	測定結果	118	135	69	81	142	0	159	0	271	285	132	136	92	0	245	0	MH025066,MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	1	1	1	1	0	1	0	1	1	1	1	1	1	1	1				13	16	81.3%
2回目	測定結果	117	135	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025066,MH025012	MH025066	0				
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
	フラグ1	1	1	1	1	1	0	1	0	1	1	1	1	1	1	1	1				14	16	87.5%	
(1)ー56	1回目	測定結果	117	123	68	81	141	0	165	0	271	0	136	0	72	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1						
2回目	測定結果	117	123	69	81	141	0	165	0	271	0	136	0	71	0	248	0	MH025362	MH025362	1				
	正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0							
	フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				15	16	93.8%	
(1)ー57	1回目	測定結果	135	0	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0				14	16	87.5%
		フラグ1	1	1	1	1	1	0	1	0	1	1	1	1	1	1	1	1						
2回目	測定結果	135	0	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1				
	正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
	フラグ1	1	1	1	1	1	0	1	0	1	1	1	1	1	1	1	1				14	16	87.5%	
(1)ー58	1回目	測定結果	117	123	69	0	128	141	159	0	267	270	136	0	72	0	209	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1						
2回目	測定結果	118	123	69	0	128	141	159	0	267	270	136	0	72	0	209	0	MH025064	MH025064	1				
	正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0							
	フラグ1	0	1	1	1	1	1	1	0	1	1	1	1	1	1	1								

(1)ー60	1回目	測定結果	117	0	64	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1				
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0							
		フラグ1	1	1	1	1	0	1	0	1	1	1	1	1	1	0	1	1				13	16	81.3%	
	2回目	測定結果	117	0	64	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1				
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0							
		フラグ1	1	1	1	1	0	1	0	1	1	1	1	1	1	0	1	1				13	16	81.3%	
		totalフラグ計	37	54	60	60	41	57	0	54	60	60	49	60	60	50	56	60		正解 品種数	100			16	81.3%
		91	120	98	54	120	109	110	116																
		一致率	38%	50%	41%	23%	50%	45%	46%	48%															

ブナシメジ品種識別検査結果票(B社・許容範囲±3)

サンプル 番号	検査	項目	msHm-1		msHm-2		msHm-3		msHm-4		msHm-5		msHm-6		msHm-7		msHm-8		判定品種	正解品種名	品種名一致	一致 フラグメント数	最大 フラグメント数	一致率
(1)ー1	1回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
(1)ー2	2回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー3	1回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	1	1	1	1	1	1	0	1	1	0	1	1	1	1	1				10	16	62.5%
(1)ー4	2回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	1	1	1	1	1	1	0	1	1	0	1	1	1	1	1				10	16	62.5%
(1)ー5	1回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1				12	16	75.0%
(1)ー6	2回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1				12	16	75.0%
(1)ー7	1回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー8	2回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー9	1回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー10	2回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー11	1回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー12	2回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー13	1回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー14	2回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー15	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー16	2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー17	1回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー18	2回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー19	1回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	1	1	1	1	1	1	0	1	1	0	1	1	1	1	1				10	16	62.5%
(1)ー20	2回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	1	1	1	1	1	1	0	1	1	0	1	1	1	1	1				10	16	62.5%
(1)ー21	1回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー22	2回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%

(1)ー12	1回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー13	2回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー14	1回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー15	2回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー16	1回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー17	2回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー18	1回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー19	2回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー20	1回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー21	2回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー22	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー23	2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー24	1回目	測定結果	117	133	66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				10	16	62.5%
(1)ー25	2回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				10	16	62.5%
(1)ー26	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー27	2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー28	1回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー29	2回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー30	1回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー31	2回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー32	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー33	2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%

(1)ー24	1回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1				12	16	75.0%
(1)ー25	2回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1				12	16	75.0%
(1)ー26	1回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー27	2回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー28	1回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー29	2回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー30	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー31	2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー32	1回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー33	2回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー34	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー35	2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー36	1回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー37	2回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー38	1回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー39	2回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				12	16	75.0%
(1)ー40	1回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%
(1)ー41	2回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				12	16	75.0%

(1)ー36	1回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1					
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0					12	16	75.0%	
		フラグ1		1		1	1	1		1		1	1	1	1		1	1								
2回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1						
	正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0									
	フラグ1		1		1	1	1		1		1	1	1	1		1	1					12	16	75.0%		
(1)ー37	1回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1					
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0					12	16	75.0%	
		フラグ1		1		1	1	1		1		1	1	1		0	1	1								
2回目	測定結果	117		62	79	145		150		270	284	131		70		243		MH025397	MH025397	1						
	正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						12	16	75.0%	
	フラグ1		1		1	1	1		1		1	1	1		0	1	1					12	16	75.0%		
(1)ー38	1回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1					
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						12	16	75.0%
		フラグ1		1		1	1	1		1		1	1	1	1		0	1					12	16	75.0%	
2回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1						
	正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							12	16	75.0%
	フラグ1		1		1	1	1		1		1	1	1	1		0	1						12	16	75.0%	
(1)ー39	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0					
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						12	16	75.0%
		フラグ1		1		1	1	1		1		1	1	1	1		1	1					12	16	75.0%	
2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025012	0						
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							12	16	75.0%
	フラグ1		1		1	1	1		1		1	1	1	1		1	1						12	16	75.0%	
(1)ー40	1回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1					
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						12	16	75.0%
		フラグ1		1		1	1	1		1		1	1	1	1		1	1						12	16	75.0%
2回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1						
	正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							12	16	75.0%
	フラグ1		1		1	1	1		1		1	1	1	1		1	1						12	16	75.0%	
(1)ー41	1回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1					
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						12	16	75.0%
		フラグ1		1		1	1	1		1		1	1	1	1		0	1					12	16	75.0%	
2回目	測定結果	123	133	66		127		150		266	284	142		90				MH025288	MH025288	1						
	正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							12	16	75.0%
	フラグ1		1		1	1	1		1		1	1	1	1		0	1						12	16	75.0%	
(1)ー42	1回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1					
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						12	16	75.0%
		フラグ1		1		1	1	1	1		1	1	1	1	1		1	1					12	16	75.0%	
2回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1						
	正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						12	16	75.0%	
	フラグ1		1		1	1	1	1		1	1	1	1	1		1	1					12	16	75.0%		
(1)ー43	1回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1					
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						12	16	75.0%
		フラグ1		1		1	1	1		1		1	1	1	1		1	1					12	16	75.0%	
2回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1						
	正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247							12	16	75.0%
	フラグ1		1		1	1	1		1		1	1	1	1		1	1						12	16	75.0%	
(1)ー44	1回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1					
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						12	16	75.0%
		フラグ1		1		1	1	1		1		1	1	1	1		1	1					12	16	75.0%	
2回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1						
	正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0							12	16	75.0%
	フラグ1		1		1	1	1		1		1	1	1	1		1	1						12	16	75.0%	
(1)ー45	1回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1					
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						12	16	75.0%
		フラグ1		1		1	1	1	1		1	1	1	1	1		1	1					12	16	75.0%	
2回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1						
	正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0							12	16	75.0%
	フラグ1		1		1	1	1	1		1	1	1	1	1		1	1						12	16	75.0%	
(1)ー46	1回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1					
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						10	16	62.5%
		フラグ1		1		1	1	1	1	0		1	0	1	1		1	1					10	16	62.5%	
2回目	測定結果	123		66	71	127	140	150	162	270		135		70		218		MH025462	MH025462	1						
	正解値	121	0	67	71	127	139	161	0	269																

(1)ー48	1回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1					
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1								
2回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1						
	正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						12	16	75.0%	
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1									
(1)ー49	1回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1					
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1								
2回目	測定結果	117	123	66	71	140		162		270	281	135		70		218	247	MH025369	MH025369	1						
	正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						12	16	75.0%	
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1									
(1)ー50	1回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1					
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1								
2回目	測定結果	133		66	79	127	140	156		270	292	131	135	70		243		MH025286	MH025286	1						
	正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0									
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1									
(1)ー51	1回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1					
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1								
2回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1						
	正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						12	16	75.0%	
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1									
(1)ー52	1回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1					
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1								
2回目	測定結果	133		66	79	145		156		270		131		70	132	243		ホクト白1号菌	ホクト白1号菌	1						
	正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0							12	16	75.0%
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1									
(1)ー53	1回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1					
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1								
2回目	測定結果	123		71	79	140		156		270	292	135		70		247		MH025463	MH025463	1						
	正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						12	16	75.0%	
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1									
(1)ー54	1回目	測定結果	123	133	66		127	140	150		266	284	142		90				MH025288	MH025288	1					
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1								
2回目	測定結果	123	133	66		127	140	150		266	284	142		90				MH025288	MH025288	1						
	正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						12	16	75.0%	
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1									
(1)ー55	1回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0					
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1								
2回目	測定結果	117	133	66	79	140		156		270	284	131	135	90		243		MH025066, MH025012	MH025066	0						
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						12	16	75.0%	
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1									
(1)ー56	1回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1					
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1								
2回目	測定結果	117	123	66	79	140		162		270		135		70		247		MH025362	MH025362	1						
	正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						12	16	75.0%	
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1									
(1)ー57	1回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1					
		正解値	133	0	67	79	140	0	155	0	269	284	130	134	92	0	243	0						6	16	37.5%
		フラグ1	1	1	0	0	1	1	1	1	0	0	0	0	0	1	1	1								
2回目	測定結果	133		66	79	140		156		270	284	131	135	90		243		MH025528	MH025528	1						
	正解値	133	0	67	79	140	0	155	0	269	284	130	134	92	0	243	0						6	16	37.5%	
	フラグ1	1	1	0	0	1	1	1	1	0	0	0	0	0	1	1	1									
(1)ー58	1回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1					
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						12	16	75.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1								
2回目	測定結果	117	123	66		127	140	156		266	270	135		70		207		MH025064	MH025064	1						
	正解値	115	121	67	0	127	139	155	0	265	269	134	0	72												

(1)ー60	1回目	測定結果	117	62	79	145	150	270	284	131	70	243	MH025397	MH025397	1				
		正解値 フラグ1	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0	
	2回目	測定結果	117	62	79	145	150	270	284	131	70	243	MH025397	MH025397	1		12	16	75.0%
		正解値 フラグ1	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0	
	totalフラグ計		60	60	60	60	60	60	60	54	60	60	54	60	60	50	56	60	
	一致率		50%	50%	50%	50%	48%	50%	48%	46%	48%								
														正解 品種数	100		平均 マーカ ー一致率	73.5%	

ブナシメジ品種識別検査結果票(F社・許容範囲±3)

サンプル 番号	検査	項目	msHm-1		msHm-2		msHm-3		msHm-4		msHm-5		msHm-6		msHm-7		msHm-8		判定品種	正解品種名	品種名一致	一致 フラグメント数	最大 フラグメント数	一致率
(1)-1	1回目	測定結果	117	122	66	71	139	0	161	0	270	281	131	135	70	0	217	247	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1	1				15	16	93.8%
(1)-1	2回目	測定結果	117	122	66	71	139	0	161	0	270	281	135	0	70	0	217	246	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
(1)-2	1回目	測定結果	122	146	66	71	127	139	149	162	270	0	131	135	70	0	217	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	0	1	1	1	1	1	0	1	1	1	1	1	1	1	1				14	16	87.5%
(1)-2	2回目	測定結果	122	146	66	71	127	139	149	161	270	0	131	135	70	0	217	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	0	1	1	1	1	1	0	1	1	1	1	1	1	1	1				14	16	87.5%
(1)-3	1回目	測定結果	122	134	66	0	127	145	149	155	266	284	131	0	70	90	217	247	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	1	1	1	1	1	0	1	0	1	1	0	1	1	0	1	0				11	16	68.8%
(1)-3	2回目	測定結果	122	134	66	0	127	0	149	0	266	284	131	143	70	90	×	0	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	0	1	0	0	1	1				13	16	81.3%
(1)-4	1回目	測定結果	117	122	66	70	139	0	162	0	270	281	131	135	70	0	217	247	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1	1				15	16	93.8%
(1)-4	2回目	測定結果	117	122	66	71	139	0	161	0	270	281	135	0	70	0	217	247	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
(1)-5	1回目	測定結果	117	122	66	79	127	139	162	0	270	0	135	0	70	0	247	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)-5	2回目	測定結果	117	122	66	79	139	0	161	0	270	0	135	0	70	0	246	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
(1)-6	1回目	測定結果	134	0	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
(1)-6	2回目	測定結果	134	0	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
(1)-7	1回目	測定結果	117	0	×	0	146	0	149	0	270	284	131	0	×	0	243	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	1	1	0	0	0	1	1	1	1	1	1	1	0	0	1	1				12	16	75.0%
(1)-7	2回目	測定結果	117	0	62	79	146	0	149	0	270	284	131	0	70	0	243	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1				15	16	93.8%
(1)-8	1回目	測定結果	117	134	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025066・MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
(1)-8	2回目	測定結果	117	134	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025066・MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
(1)-9	1回目	測定結果	134	146	66	79	127	139	156	0	270	292	131	0	70	0	243	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)-9	2回目	測定結果	134	146	67	79	127	139	155	0	270	292	131	0	70	0	242	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)-10	1回目	測定結果	122	146	66	71	127	139	149	161	270	0	131	135	70	0	217	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	0	1	1	1	1	1	0	1	1	1	1	1	1	1	1				14	16	87.5%
(1)-10	2回目	測定結果	117	122	66	71	127	139	149	161	270	0	131	135	70	0	217	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	0	1	1	1	1	1	0	1	1	1	1	1	1	1	1				14	16	87.5%
(1)-11	1回目	測定結果	117	122	66	0	127	139	156	0	266	270	131	135	70	0	206	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	0	1	1	1	1	1				15	16	93.8%
(1)-11	2回目	測定結果	117	122	66	0	127	139	156	0	266	270	135	143	70	0	×	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	0	1	1	0	1	1				14	16	87.5%

(1)ー12	1回目	測定結果	134	0	66	79	145	0	156	0	270	0	131	0	70	132	243	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				15	16	93.8%
2回目	測定結果	134	0	66	79	145	0	155	0	270	0	131	0	×	0	243	0	ホクト白1号菌	ホクト白1号菌	1				
	正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0							
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1				15	16	93.8%	
(1)ー13	1回目	測定結果	116	122	66	79	139	0	161	0	270	0	135	0	70	0	246	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
2回目	測定結果	117	122	66	79	139	0	161	0	270	0	135	0	70	0	246	0	MH025362	MH025362	1				
	正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0							
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%	
(1)ー14	1回目	測定結果	122	146	71	79	139	0	156	0	270	292	131	135	70	0	247	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	0	1	1	1	1				14	16	87.5%
2回目	測定結果	122	146	71	79	139	0	155	0	×	0	135	0	70	0	247	0	MH025463	MH025463	1				
	正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0							
	フラグ1	1	0	1	1	1	1	1	1	0	0	1	1	1	1	1	1				13	16	81.3%	
(1)ー15	1回目	測定結果	134	0	66	79	145	0	155	0	270	0	131	0	70	132	243	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				15	16	93.8%
2回目	測定結果	134	0	66	79	145	0	155	0	270	0	131	0	70	132	242	0	ホクト白1号菌	ホクト白1号菌	1				
	正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0							
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				15	16	93.8%	
(1)ー16	1回目	測定結果	134	0	66	79	127	139	155	0	270	292	131	0	70	0	243	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
2回目	測定結果	134	146	67	79	127	139	156	0	270	292	131	135	70	0	243	0	MH025286	MH025286	1				
	正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0							
	フラグ1	1	0	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%	
(1)ー17	1回目	測定結果	×	0	66	79	139	0	155	0	270	284	131	135	90	0	243	0	MH025066・MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	0	0	1	1	1	1	1	1	1	1	1	1	1	1	1	1				14	16	87.5%
2回目	測定結果	117	134	66	79	139	0	155	0	270	284	131	135	90	0	243	0	MH025066・MH025012	MH025066	0				
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%	
(1)ー18	1回目	測定結果	122	146	66	71	127	139	149	161	270	0	131	135	70	0	217	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	0	1	1	1	1	1	0	1	1	1	1	1	1	1	1				14	16	87.5%
2回目	測定結果	122	146	66	71	127	139	149	161	270	0	131	135	70	0	217	0	MH025462	MH025462	1				
	正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0							
	フラグ1	1	0	1	1	1	1	1	0	1	1	1	1	1	1	1	1				14	16	87.5%	
(1)ー19	1回目	測定結果	117	134	66	79	139	0	155	0	270	284	131	135	90	0	242	0	MH025066・MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
2回目	測定結果	117	134	66	79	139	0	156	0	270	284	131	135	90	0	242	0	MH025066・MH025012	MH025012	0				
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%	
(1)ー20	1回目	測定結果	117	134	66	79	139	0	155	0	270	284	131	135	70	90	246	0	MH025066・MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	0	1	1	1	1	1	0	1	0	0	1	0	1	1	1				11	16	68.8%
2回目	測定結果	117	134	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025066・MH025012	MH025012	0				
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
	フラグ1	1	0	1	1	1	1	1	0	1	0	0	1	0	1	0	1				10	16	62.5%	
(1)ー21	1回目	測定結果	117	122	66	71	139	145	156	161	270	281	131	0	70	0	217	246	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	0	1	0	1	0	0	1	1	1	1	0	1	1	0	0				9	16	56.3%
2回目	測定結果	117	122	66	71	139	0	161	0	270	281	135	0	70	132	217	247	MH025369	MH025369	1				
	正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247							
	フラグ1	1	0	1	0	1	1	0	1	1	1	1	0	0	1	0	0				9	16	56.3%	
(1)ー22	1回目	測定結果	134	0	66	79	145	0	156	0	270	0	131	0	70	132	243	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	0	1	0	0	1	1	1				12	16	75.0%
2回目	測定結果	134	0	66	79	145	0	155	0	270	0	131	0	70	132	243	0	ホクト白1号菌	ホクト白1号菌	1				
	正解値	133	0	67	79	145	0	155																

(1)ー24	1回目	測定結果	×	0	66	0	127	0	149	156	266	284	131	0	70	90	221	0	MH025288	MH025288	1				
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
		フラグ1	0	0	1	1	1	1	1	0	1	1	0	1	1	0	1	1					11	16	68.8%
2回目	測定結果	122	134	66	0	127	0	149	156	266	284	131	143	70	90	×	0	MH025288	MH025288	1					
	正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0								
	フラグ1	1	1	1	1	1	1	1	0	1	1	1	0	1	0	0	1					12	16	75.0%	
(1)ー25	1回目	測定結果	117	122	66	0	127	139	155	0	266	270	131	135	70	0	246	0	MH025064	MH025064	1				
		正解値	115	121	67	0	127	139	155	0	265	269	134	72	0	207	0								
		フラグ1	1	0	1	0	0	1	1	1	0	1	1	0	0	1	0	1					9	16	56.3%
2回目	測定結果	117	122	66	0	127	139	155	0	×	0	135	0	70	0	206	0	MH025064	MH025064	1					
	正解値	115	121	67	0	127	139	155	0	265	269	134	72	0	207	0									
	フラグ1	1	0	1	0	0	1	1	1	0	1	1	0	0	1	0	1					9	16	56.3%	
(1)ー26	1回目	測定結果	134	0	66	79	139	0	155	0	270	284	131	135	90	0	243	0	MH025528	MH025528	1				
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0							
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1					15	16	93.8%
2回目	測定結果	134	0	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025528	MH025528	1					
	正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0								
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1					15	16	93.8%	
(1)ー27	1回目	測定結果	117	134	66	79	139	145	155	0	270	284	131	135	70	0	242	0	MH025066・MH025012	MH025066	0				
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	1	0	1	0	1	0	1	0	1	1	1	1	0	1	1	0					10	16	62.5%
2回目	測定結果	117	134	66	79	139	0	156	0	270	284	131	135	90	0	242	0	MH025066・MH025012	MH025066	0					
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0								
	フラグ1	1	0	1	0	1	1	0	1	0	1	0	1	0	0	1	0					8	16	50.0%	
(1)ー28	1回目	測定結果	116	122	66	79	139	0	161	0	270	0	135	0	70	0	246	0	MH025362	MH025362	1				
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0							
		フラグ1	1	0	1	1	1	1	0	1	1	0	1	0	1	0	1	1					11	16	68.8%
2回目	測定結果	116	122	66	79	139	0	161	0	270	0	135	0	70	0	247	0	MH025362	MH025362	1					
	正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0								
	フラグ1	1	0	1	1	1	1	0	1	0	1	0	0	1	0	1	0					10	16	62.5%	
(1)ー29	1回目	測定結果	117	134	66	79	139	145	155	0	270	284	131	135	70	90	243	0	MH025066・MH025012	MH025012	0				
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	1	1	1	1	1	0	1	1	1	1	1	1	1	0	1	1					14	16	87.5%
2回目	測定結果	117	134	66	79	139	0	155	0	270	284	131	135	90	0	242	0	MH025066・MH025012	MH025012	0					
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0								

(1)ー36	1回目	測定結果	134	0	66	79	145	0	156	0	270	284	131	135	90	0	242	0	MH025528	MH025528	1				
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
	2回目	測定結果	134	0	66	79	139	0	156	0	270	284	131	135	90	0	242	0	MH025528	MH025528	1				
(1)ー37	1回目	正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
		測定結果	117	0	62	79	146	0	149	0	270	284	131	0	70	0	242	0	MH025397	MH025397	1				
	2回目	正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0							
(1)ー38	1回目	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
		測定結果	116	0	62	79	146	0	149	0	270	284	131	0	70	0	243	0	MH025397	MH025397	1				
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0							
	2回目	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1	1				15	16	93.8%
(1)ー39	1回目	測定結果	117	122	66	0	127	0	155	0	266	284	131	0	90	0	×	0	MH025288	MH025288	1				
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
		フラグ1	1	0	1	1	1	1	1	0	1	1	0	1	1	1	0	1	1				12	16	75.0%
	2回目	測定結果	122	134	66	0	127	0	149	0	266	284	131	143	70	90	243	0	MH025288	MH025288	1				
(1)ー40	1回目	正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
		フラグ1	1	1	1	1	1	1	1	1	1	1	0	1	0	0	1	1	1				13	16	81.3%
		測定結果	111	117	66	79	139	145	156	0	270	284	131	135	90	0	242	0	MH025066・MH025012	MH025012	0				
	2回目	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
(1)ー41	1回目	フラグ1	1	0	1	1	1	0	1	1	1	1	1	1	1	1	1	1	1				14	16	87.5%
		測定結果	117	134	66	79	139	0	155	0	270	284	131	135	90	0	243	0	MH025066・MH025012	MH025012	0				
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
	2回目	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
(1)ー42	1回目	測定結果	134	0	66	79	139	0	156	0	270	284	135	0	70	90	242	0	MH025528	MH025528	1				
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	1	1	1	1	1	1	1	1	1	1	0	1	1	0	1	1	1				14	16	87.5%
	2回目	測定結果	134	0	66	79	139	0	156	0	270	284	131	135	90	0	243	0	MH025528	MH025528	1				
(1)ー43	1回目	正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				16	16	100.0%
		測定結果	122	134	66	0	128	0	149	155	266	284	131	0	70	90	206	221	MH025288	MH025288	1				
	2回目	正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
(1)ー44	1回目	フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	0	1	1	1				12	16	75.0%
		測定結果	122	134	66	0	127	0	149	155	266	284	143	0	70	90	×	0	MH025288	MH025288	1				
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
	2回目	フラグ1	1	1	1	1	1	1	1	0	1	1	1	1	1	0	0	1	1				13	16	81.3%
(1)ー45	1回目	測定結果	134	0	66	79	127	139	156	0	270	292	131	0											

(1)－48	1回目	測定結果	122	146	71	79	139	0	156	0	270	292	135	0	70	0	246	0	MH025463	MH025463	1				
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0					15	16	93.8%
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	1	1	1	1	1							
		測定結果	122	146	71	79	139	0	155	0	270	292	135	0	70	0	246	0	MH025463	MH025463	1				
2回目	1回目	正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0							
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	1	1	1	1	1					15	16	93.8%
		測定結果	117	122	66	71	139	0	161	0	270	281	135	0	70	0	217	246	MH025369	MH025369	1				
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247					16	16	100.0%
(1)－49	2回目	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1					16	16	100.0%
		測定結果	117	122	66	71	139	0	161	0	270	281	135	0	70	0	217	246	MH025369	MH025369	1				
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247					16	16	100.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1					16	16
(1)－50	1回目	測定結果	134	0	66	79	127	139	155	0	270	292	131	135	70	0	242	0	MH025286	MH025286	1				
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0					16	16	100.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1					16	16	100.0%
		測定結果	134	146	67	79	127	139	156	0	270	292	131	0	70	0	217	242	MH025286	MH025286	1				
(1)－51	2回目	正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0							
		フラグ1	1	0	1	1	1	1	1	1	1	1	1	1	1	1	1	0					14	16	87.5%
		測定結果	117	122	66	79	139	0	161	0	270	0	135	0	70	0	246	0	MH025362	MH025362	1				
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0					16	16	100.0%
(1)－52	2回目	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1					16	16	100.0%
		測定結果	116	122	66	79	139	0	161	0	270	0	135	0	70	0	246	0	MH025362	MH025362	1				
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0					16	16	100.0%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1					16	16
(1)－53	1回目	測定結果	122	134	66	79	145	0	156	0	270	0	131	135	70	132	242	0	ホクト白1号菌	ホクト白1号菌	1				
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0							
		フラグ1	1	0	1	1	1	1	1	1	1	1	0	1	0	1	0	1					13	16	81.3%
		測定結果	122	134	66	79	145	0	156	0	270	0	131	135	70	132	242	246	ホクト白1号菌	ホクト白1号菌	1				
(1)－54	2回目	正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0							
		フラグ1	1	0	1	1	1	1	1	1	1	1	0	1	0	1	0	1					12	16	75.0%
		測定結果	122	146	71	79	139	0	156	0	270	292	135	0	70	0	246	0	MH025463	MH025463	1				
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0					15	16	93.8%
(1)－55	2回目	フラグ1	1	0	1	1	1	1	1	1	1	1	1	1	1	1	1	1					15	16	93.8%
		測定結果	122	146	71	79	127	139	155	0	270	292	135	0	70	132	246	0	MH025463	MH025463					

(1)～60	1回目	測定結果	117	0	62	79	127	146	149	0	270	284	131	0	70	0	242	0	MH025397	MH025397	1				
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0							
		フラグ1	1	1	1	1	1	0	1	1	1	1	1	1	1	0	1	1					14	16	87.5%
	2回目	測定結果	117	0	62	79	146	0	149	0	270	284	131	0	70	0	242	0	MH025397	MH025397	1				
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0							
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1					15	16	93.8%
	totalフラグ計	58	35	59	51	56	52	53	48	56	51	54	45	46	44	49	55		正解 品種数	100			平均 マーカ ー致率	88.1%	
		93		110		108		101		107		99		90		104									
	一致率	39%		46%		45%		42%		45%		41%		38%		43%									

ブナシメジ品種識別検査結果票(L社・許容範囲±3)

サンプル番号	検査	項目	msHm-1		msHm-2		msHm-3		msHm-4		msHm-5		msHm-6		msHm-7		msHm-8		判定品種	正解品種名	品種名一致	一致フラグメント数	最大フラグメント数	一致率
(1)-1	1回目	測定結果	117	123	69	73	141	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				15	16	93.8%
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1						
(1)-2	2回目	測定結果	118	124	69	73	142	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				15	16	93.8%
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1						
(1)-3	1回目	測定結果	123	0	69	73	128	142	152	165	270	0	136	0	72	0	219	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0				13	16	81.3%
		フラグ1	1	1	1	1	1	1	0	0	1	1	1	1	1	1	1	1						
(1)-4	2回目	測定結果	124	0	69	73	128	141	152	165	271	0	136	0	72	0	219	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0				13	16	81.3%
		フラグ1	1	1	1	1	1	1	0	0	1	1	1	1	1	1	1	1						
(1)-5	1回目	測定結果	123	135	69	0	128	0	152	0	267	285	145	0	92	0	0	0	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1						
(1)-6	2回目	測定結果	123	135	69	0	128	0	152	0	267	285	146	0	92	0	0	0	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0				14	16	87.5%
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1						
(1)-7	1回目	測定結果	118	124	69	73	141	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				15	16	93.8%
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1						
(1)-8	2回目	測定結果	118	124	69	73	141	0	165	0	271	282	136	0	71	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247				15	16	93.8%
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1						
(1)-9	1回目	測定結果	117	124	69	81	141	0	165	0	271	0	137	0	72	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1						
(1)-10	2回目	測定結果	118	124	69	81	141	0	165	0	271	0	136	0	71	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1						
(1)-11	1回目	測定結果	135	0	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1						
(1)-12	2回目	測定結果	135	0	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1						
(1)-13	1回目	測定結果	118	0	65	81	149	0	152	0	271	285	132	0	72	0	245	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0				14	16	87.5%
		フラグ1	1	1	1	1	0	1	1	1	1	1	1	1	1	0	1	1						
(1)-14	2回目	測定結果	118	0	65	81	149	0	153	0	271	285	132	0	71	0	244	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0				13	16	81.3%
		フラグ1	1	1	1	1	0	1	0	1	1	1	1	1	1	0	1	1						
(1)-15	1回目	測定結果	118	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1						
(1)-16	2回目	測定結果	118	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1						
(1)-17	1回目	測定結果	135	0	69	81	128	141	159	0	271	293	132	0	72	0	244	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1						
(1)-18	2回目	測定結果	135	0	69	81	129	141	159	0	271	293	132	0	72	0	244	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1						
(1)-19	1回目	測定結果	123	0	69	73	128	141	152	165	271	0	136	0	72	0	219	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0				13	16	81.3%
		フラグ1	1	1	1	1	1	1	0	0	1	1	1	1	1	1	1	1						
(1)-20	2回目	測定結果	123	0	69	73	128	142	152	165	270	0	136	0	71	0	219	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0				13	16	81.3%
		フラグ1	1	1	1	1	1	1	0	0	1	1	1	1	1	1	1	1						
(1)-21	1回目	測定結果	117	123	69	0	128	141	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1						
(1)-22	2回目	測定結果	118	123	69	0	128	142	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0				15	16	93.8%
		フラグ1	1	1	1	1	1	1	0	1	1	1	1	1	1	1	1	1						

(1)ー12	1回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	72	133	244	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				14	16	87.5%
(1)ー13	2回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	72	133	244	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				14	16	87.5%
(1)ー14	1回目	測定結果	117	124	69	81	141	0	165	0	271	0	136	0	72	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー15	2回目	測定結果	118	123	69	81	141	0	165	0	271	0	136	0	71	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー16	1回目	測定結果	123	0	73	81	141	0	159	0	271	293	136	0	72	0	248	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー17	2回目	測定結果	124	0	73	81	141	0	159	0	271	293	136	0	71	0	248	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー18	1回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	72	133	244	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				14	16	87.5%
(1)ー19	2回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	72	133	244	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				14	16	87.5%
(1)ー20	1回目	測定結果	135	0	69	81	128	141	159	0	271	293	132	0	72	0	244	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー21	2回目	測定結果	135	0	69	81	129	141	159	0	271	293	132	0	72	0	244	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー22	1回目	測定結果	117	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066,MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー23	2回目	測定結果	118	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066,MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー24	1回目	測定結果	123	0	69	72	128	141	152	165	271	0	136	0	72	0	219	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				13	16	81.3%
(1)ー25	2回目	測定結果	123	0	69	72	128	141	152	165	271	0	136	0	71	0	219	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				13	16	81.3%
(1)ー26	1回目	測定結果	117	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066,MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー27	2回目	測定結果	118	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066,MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー28	1回目	測定結果	118	123	69	73	142	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー29	2回目	測定結果	118	123	69	73	141	0	165	0	271	282	136	0	71	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー30	1回目	測定結果	135	0	68	81	148	0	159	0	271	0	132	0	72	133	244	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				14	16	87.5%
(1)ー31	2回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	71	134	244	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				14	16	87.5%
(1)ー32	1回目	測定結果	117	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066,MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー33	2回目	測定結果	118	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066,MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%

(1)ー24	1回目	測定結果	123	135	69	0	128	0	152	0	267	285	145	0	92	0	244	0	MH025288	MH025288	1					
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0								
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1				15	16	93.8%	
2回目	測定結果	124	135	68	0	128	0	152	0	267	285	145	0	92	0	244	0	MH025288	MH025288	1						
	正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0									
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1				15	16	93.8%		
(1)ー25	1回目	測定結果	117	123	69	0	128	141	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1					
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0								
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%		
2回目	測定結果	117	123	68	0	128	141	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1						
	正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0									
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%			
(1)ー26	1回目	測定結果	135	0	68	81	141	0	159	0	271	285	132	136	92	0	245	0	MH025528	MH025528	1					
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0								
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%		
2回目	測定結果	135	0	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1						
	正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0									
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%			
(1)ー27	1回目	測定結果	117	135	68	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025066	0					
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0								
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%		
2回目	測定結果	117	135	69	81	142	0	159	0	270	285	132	136	92	0	244	0	MH025066.MH025012	MH025066	0						
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0									
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%			
(1)ー28	1回目	測定結果	117	123	68	81	142	0	165	0	271	0	136	0	72	0	248	0	MH025362	MH025362	1					
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0								
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%		
2回目	測定結果	118	123	69	81	141	0	165	0	271	0	136	0	71	0	248	0	MH025362	MH025362	1						
	正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0									
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%			
(1)ー29	1回目	測定結果	117	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025012	0					
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0								
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%		
2回目	測定結果	118	135	69	81	142	0	159	0	271	285	132	137	92	0	244	0	MH025066.MH025012	MH025012	0						
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0									
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%			
(1)ー30	1回目	測定結果	118	0	64	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1					
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0								
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1				14	16	87.5%		
2回目	測定結果	118	0	65	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1						
	正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0									
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1				14	16	87.5%			
(1)ー31	1回目	測定結果	117	123	69	0	128	141	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1					
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0								
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%		
2回目	測定結果	117	123	69	0	128	141	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1						
	正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0									
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%			
(1)ー32	1回目	測定結果	123	0	72	81	141	0	159	0	271	293	136	0	72	0	248	0	MH025463	MH025463	1					
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0								
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%		
2回目	測定結果	123	0	73	81	141	0	159	0	271	293	136	0	71	0	248	0	MH025463	MH025463	1						
	正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0									
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%			
(1)ー33	1回目	測定結果	117	0	65	81	149	0	152	0	271	285	132	0	72	0	245	0	MH025397	MH025397	1					
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0								
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				14	16	87.5%		
2回目	測定結果	117	0	65	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1						
	正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0									
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				14	16	87.5%			
(1)ー34	1回目	測定結果	117	123	69	0	128	141	159	0	267	271	136	0	72	0	209	0	MH025064	MH025064	1					
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0								
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%		
2回目	測定結果	117	123	69																						

(1)ー36	1回目	測定結果	135	0	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー37	2回目	測定結果	135	0	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー38	1回目	測定結果	117	0	64	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	1	1	1	1	0	1	1	1	1	1	1	1	1	0	1	1				14	16	87.5%
(1)ー39	2回目	測定結果	117	0	65	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1			
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0						
		フラグ1	1	1	1	1	0	1	1	1	1	1	1	1	1	0	1	1				14	16	87.5%
(1)ー40	1回目	測定結果	123	135	69	0	128	0	152	0	267	285	145	0	92	0	245	0	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1				15	16	93.8%
(1)ー41	2回目	測定結果	123	135	69	0	128	0	152	0	267	285	145	0	92	0	244	0	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1				15	16	93.8%
(1)ー42	1回目	測定結果	117	135	69	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー43	2回目	測定結果	117	135	69	81	141	0	159	0	271	285	132	136	92	0	245	0	MH025066.MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー44	1回目	測定結果	135	0	69	81	141	0	159	0	271	285	132	136	92	0	245	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー45	2回目	測定結果	135	0	69	81	141	0	159	0	271	285	132	136	92	0	245	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー46	1回目	測定結果	135	0	69	81	128	141	159	0	271	293	132	0	72	0	244	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー47	2回目	測定結果	135	0	69	81	128	142	159	0	270	293	132	0	72	0	244	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー48	1回目	測定結果	117	123	68	73	142	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー49	2回目	測定結果	118	124	69	73	141	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー50	1回目	測定結果	123	0	73	81	141	0	159	0	271	293	136	0	72	0	248	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー51	2回目	測定結果	124	0	73	81	141	0	159	0	271	293	136	0	71	0	248	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー52	1回目	測定結果	135	0	69	81	128	141	159	0	271	293	132	0	72	0	244	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー53	2回目	測定結果	135	0	69	81	128	141	159	0	271	293	132	0	71	0	244	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー54	1回目	測定結果	124	0	69	73	128	141	152	165	271	0	136	0	72	0	219	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	1	1	1	1	1	1	0	0	1	1	1	1	1	1	1				13	16	81.3%
(1)ー55	2回目	測定結果	124	0	69	73	128	141	152	165	271	0	136	0	72	0	219	0	MH025462	MH025462	1			
		正解値	121	0	67	71	127	139	161	0	269	0	130	134	72	0	217	0						
		フラグ1	1	1	1	1	1	1	1	0	0	1	1	1	1	1	1	1				13	16	81.3%
(1)ー56	1回目	測定結果	118	135	68	81	141	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
(1)ー57	2回目	測定結果	118	135	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025012	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%

(1)ー48	1回目	測定結果	123	0	72	81	141	0	159	0	271	293	136	0	72	0	248	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
2回目	測定結果	123	0	72	81	142	0	159	0	271	293	136	0	71	0	248	0	MH025463	MH025463	1				
	正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0							
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%	
(1)ー49	1回目	測定結果	118	124	68	72	141	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1			
		正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
2回目	測定結果	118	123	69	72	141	0	165	0	271	282	136	0	72	0	219	248	MH025369	MH025369	1				
	正解値	115	121	67	71	139	0	161	0	269	281	134	0	72	0	217	247							
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%	
(1)ー50	1回目	測定結果	135	0	69	81	128	141	159	0	271	293	132	0	72	0	244	0	MH025286	MH025286	1			
		正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
2回目	測定結果	135	0	69	81	128	141	159	0	271	292	132	0	72	0	245	0	MH025286	MH025286	1				
	正解値	133	0	67	79	127	139	155	0	269	293	130	134	72	0	243	0							
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%	
(1)ー51	1回目	測定結果	117	123	68	81	141	0	165	0	271	0	136	0	72	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
2回目	測定結果	118	123	69	81	142	0	165	0	271	0	136	0	72	0	248	0	MH025362	MH025362	1				
	正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0							
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%	
(1)ー52	1回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	72	133	244	0	ホクト白1号菌	ホクト白1号菌	1			
		正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				14	16	87.5%
2回目	測定結果	135	0	69	81	148	0	159	0	271	0	132	0	72	134	244	0	ホクト白1号菌	ホクト白1号菌	1				
	正解値	133	0	67	79	145	0	155	0	269	0	130	0	132	0	243	0							
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1	1				14	16	87.5%	
(1)ー53	1回目	測定結果	124	0	72	81	141	0	159	0	271	293	136	0	72	0	248	0	MH025463	MH025463	1			
		正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
2回目	測定結果	124	0	73	81	141	0	159	0	271	293	136	0	72	0	248	0	MH025463	MH025463	1				
	正解値	121	0	71	79	139	0	155	0	269	293	134	0	72	0	247	0							
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%	
(1)ー54	1回目	測定結果	123	135	69	0	128	0	152	0	267	285	145	0	92	0	248	0	MH025288	MH025288	1			
		正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1				15	16	93.8%
2回目	測定結果	123	135	69	0	128	0	152	0	267	285	145	0	92	0	248	0	MH025288	MH025288	1				
	正解値	121	133	67	0	127	0	149	0	265	284	142	0	92	0	220	0							
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	0	1				15	16	93.8%	
(1)ー55	1回目	測定結果	118	135	69	81	142	0	159	0	271	285	132	136	92	0	245	0	MH025066.MH025012	MH025066	0			
		正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
2回目	測定結果	117	135	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025066.MH025012	MH025066	0				
	正解値	115	133	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%	
(1)ー56	1回目	測定結果	117	123	68	81	141	0	165	0	271	0	136	0	72	0	248	0	MH025362	MH025362	1			
		正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
2回目	測定結果	117	123	69	81	141	0	165	0	271	0	136	0	71	0	248	0	MH025362	MH025362	1				
	正解値	115	121	67	79	139	0	161	0	269	0	134	0	72	0	247	0							
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%	
(1)ー57	1回目	測定結果	135	0	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1			
		正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
2回目	測定結果	135	0	69	81	142	0	159	0	271	285	132	136	92	0	244	0	MH025528	MH025528	1				
	正解値	133	0	67	79	139	0	155	0	269	284	130	134	92	0	243	0							
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%	
(1)ー58	1回目	測定結果	117	123	69	0	128	141	159	0	267	270	136	0	72	0	209	0	MH025064	MH025064	1			
		正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0						
		フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				15	16	93.8%
2回目	測定結果	118	123	69	0	128	141	159	0	267	270	136	0	72	0	209	0	MH025064	MH025064	1				
	正解値	115	121	67	0	127	139	155	0	265	269	134	0	72	0	207	0							
	フラグ1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1								

(1)ー60	1回目	測定結果	117	0	64	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1				
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0							
		フラグ1	1	1	1	1	1	0	1	1	1	1	1	1	1	0	1	1				14		16	87.5%
	2回目	測定結果	117	0	64	81	149	0	152	0	271	285	132	0	72	0	244	0	MH025397	MH025397	1				
		正解値	115	0	63	79	145	0	149	0	269	284	130	0	72	132	243	0							
		フラグ1	1	1	1	1	1	0	1	1	1	1	1	1	1	0	1	1				14		16	87.5%
		totalフラグ計	60	60	60	60	56	60	7	54	60	60	53	60	60	50	56	60		正解 品種数	100			平均 マーカ ー致率	91.6%
			一致率	50%		50%		48%		25%		50%		47%		46%		48%							
			120		120		116		61		120		113		110		116								