

添付資料 5

品種毎のマーカーー一致率

品種毎のマーカーー致率

正解品種名	最大 フラグメント 数	許容範囲	フラグメント一致数					フラグメント一致率				
			B社	F社	L社	検査毎 平均	サンプル毎 平均	B社	F社	L社	検査毎 平均	サンプル毎 平均
ホクト白1号苗	16	±1	12.0	13.8	7.9	11.2	11.2	75.0%	86.3%	49.4%	70.2%	70.2%
		±2	12.0	13.8	13.0	12.9	12.9	75.0%	86.3%	81.3%	80.8%	80.8%
		±3	12.0	13.8	14.0	13.3	13.3	75.0%	86.3%	87.5%	82.9%	82.9%
MH025066	16	±1	11.0	11.0	7.1	9.7	9.7	68.8%	68.8%	44.4%	60.6%	60.6%
		±2	12.0	12.5	14.1	12.9	12.9	75.0%	78.1%	88.1%	80.4%	80.4%
		±3	12.0	13.0	15.0	13.3	13.3	75.0%	81.3%	93.8%	83.3%	83.3%
MH025012	16	±1	11.0	12.8	7.0	10.3	10.3	68.8%	80.0%	43.8%	64.2%	64.2%
		±2	12.0	14.4	14.2	13.5	13.5	75.0%	90.0%	88.8%	84.6%	84.6%
		±3	12.0	14.5	15.0	13.8	13.8	75.0%	90.6%	93.8%	86.5%	86.5%
MH025064	16	±1	10.0	12.0	7.3	9.8	9.8	62.5%	75.0%	45.6%	61.0%	61.0%
		±2	12.0	13.7	14.7	13.5	13.5	75.0%	85.6%	91.9%	84.2%	84.2%
		±3	12.0	13.7	15.0	13.6	13.6	75.0%	85.6%	93.8%	84.8%	84.8%
MH025286	16	±1	12.0	13.6	7.8	11.1	11.1	75.0%	85.0%	48.8%	69.6%	69.6%
		±2	12.0	14.6	14.9	13.8	13.8	75.0%	91.3%	93.1%	86.5%	86.5%
		±3	12.0	15.3	15.0	14.1	14.1	75.0%	95.6%	93.8%	88.1%	88.1%
MH025288	16	±1	11.0	10.8	8.9	10.2	10.2	68.8%	67.5%	55.6%	64.0%	64.0%
		±2	12.0	11.8	12.9	12.2	12.2	75.0%	73.8%	80.6%	76.5%	76.5%
		±3	12.0	12.1	14.9	13.0	13.0	75.0%	75.6%	93.1%	81.3%	81.3%
MH025362	16	±1	10.0	13.0	8.3	10.4	10.4	62.5%	81.3%	51.9%	65.2%	65.2%
		±2	12.0	14.7	14.0	13.6	13.6	75.0%	91.9%	87.5%	84.8%	84.8%
		±3	12.0	14.7	15.0	13.9	13.9	75.0%	91.9%	93.8%	86.9%	86.9%
MH025369	16	±1	10.3	12.1	7.4	9.9	9.9	64.4%	75.6%	46.3%	62.1%	62.1%
		±2	12.4	14.0	13.4	13.3	13.3	77.5%	87.5%	83.8%	82.9%	82.9%
		±3	12.4	14.3	15.0	13.9	13.9	77.5%	89.4%	93.8%	86.9%	86.9%
MH025397	16	±1	11.0	12.5	8.2	10.6	10.6	68.8%	78.1%	51.3%	66.0%	66.0%
		±2	12.0	14.2	12.6	12.9	12.9	75.0%	88.8%	78.8%	80.8%	80.8%
		±3	12.0	14.2	13.9	13.4	13.4	75.0%	88.8%	86.9%	83.5%	83.5%
MH025462	16	±1	9.0	13.0	6.5	9.5	9.5	56.3%	81.3%	40.6%	59.4%	59.4%
		±2	10.0	14.0	12.5	12.2	12.2	62.5%	87.5%	78.1%	76.0%	76.0%
		±3	10.0	14.0	13.0	12.3	12.3	62.5%	87.5%	81.3%	77.1%	77.1%
MH025463	16	±1	11.0	13.2	9.4	11.2	11.2	68.8%	82.5%	58.8%	70.0%	70.0%
		±2	12.0	14.2	14.5	13.6	13.6	75.0%	88.8%	90.6%	84.8%	84.8%
		±3	12.0	14.2	15.0	13.7	13.7	75.0%	88.8%	93.8%	85.8%	85.8%
MH025528	16	±1	10.8	14.5	7.8	11.0	11.0	67.5%	90.6%	48.8%	69.0%	69.0%
		±2	10.8	15.3	14.5	13.5	13.5	67.5%	95.6%	90.6%	84.6%	84.6%
		±3	10.8	15.3	15.0	13.7	13.7	67.5%	95.6%	93.8%	85.6%	85.6%
平均	16	±1	10.8	12.7	7.8	10.4	10.4	67.2%	79.3%	48.8%	65.1%	65.1%
		±2	11.8	13.9	13.8	13.2	13.2	73.5%	87.1%	86.1%	82.2%	82.2%
		±3	16.0	11.8	14.1	14.7	13.5	73.5%	88.1%	91.6%	84.4%	84.4%

